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WARNING!

For a full understanding of the performance characteristics of this medical device, the user should carefully read this Instruction for Use before use of the medical device.

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For Your Safety and that of Your Patients

Strictly follow the Instructions for Use

Any use of the anaesthetic gas scavenging system requires full understanding and strict observation of these instructions.

The anaesthetic gas scavenging system is only to be used for purposes specified here.

Maintenance

The anaesthetic gas scavenging system must be inspected and serviced by trained service personnel at regular 6 month intervals.

Repair and general overhaul of the anaesthetic gas scavenging system may only be carried out by trained service personnel.

We recommend that a service contract be obtained with DrägerService and that all repairs also be carried out by them.

Only authentic Dräger spare parts may be used for maintenance.

Observe chapter "Maintenance".

Liability for proper function or damage

The liability for the proper function of the anaesthetic gas scavenging system is irrevocably transferred to the owner or operator to the extent that the anaesthetic gas scavenging system is serviced or repaired by personnel not employed or authorized by DrägerService or if the operating light is used in a manner not conforming to its intended use.

Dräger cannot be held responsible for damage caused by non-compliance with the recommendations given above.

The warranty and liability provisions of the terms of sale and delivery of Dräger are likewise not modified by the recommendations given above.

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Application

Intended Use

Anaesthetic Gas Scavenging System – AGSS

For scavenging expiratory and excess gas from anaesthetic workstations.

For reducing the concentration of anaesthetic agent in the vicinity of anaesthetic workstations.

The system is designed for the following anesthetic gases and vapours: nitrous oxide (N₂O), xenon, enflurane, isoflurane, sevoflurane, desflurane and halothane.

Description

The Anaesthetic gas scavenging system is powered by compressed air from the piped medical gas system.

When the probe is connected, it activates the supply of compressed air to ejector. The indicator changes to green to show that the anaesthetic gas scavenging system is in operation.

The surplus anaesthetic gas is exhausted from the waste gas connection of the receiving system and carried to open air or into the exhaust ducts of a non-return ventilation system.

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
Operation

Insertion

- Attach sleeve of suction hose to the waste gas connection on the receiving system.
- 1 Insert probe into terminal unit and push until it engages.
 - 2 The indicator will change to green; suction will commence.

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
Disconnection

- 3 Push release bushing and at the same time pull probe out of terminal unit.
- The probe will be released; suction will be switched off.
- 2 The indicator point will cease to be green.

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
Symboly

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Cleaning and Disinfecting

Do not use cleaning agents which release chlorine or oxygen.

- Wipe the surfaces of the terminal unit and probe with a disposable cloth.

Wipe disinfection only.

- Follow disinfectant manufacturer's instructions.
- Do not let any liquid get into the terminal unit or probe.
- Use surface disinfectants. For material compatibility use disinfectants based on
 - aldehydes,
 - alcohols,
 - quaternary ammonium compounds.

Damage to materials may occur if disinfectants are used which are based on

- phenols,
- halogen-releasing compounds,
- strong organic acids,
- oxygen-releasing compounds.

For users in Germany we recommend that only disinfectants on the current DGHM list are used (DGHM: German Society for Hygiene and Microbiology).

The DGHM list indicates the composition of each disinfectant.

For countries where the DGHM list is not available, we recommend the types of disinfectant given above.

Examples of suitable disinfectants:

- Terralin (Schülke & Mayr, Norderstedt)
- Sagromed Spray (Schülke & Mayr, Norderstedt)
- Incidin Perfect (Henkel KG aA)
- Sokrena (Bode Chemie GmbH & Co.)

Maintenance

Inspection and Maintenance	to be carried out every six month by trained service personnel
Replace sealing rings of the terminal unit	every 6 years by trained service personnel
Repairs	by trained service personnel only

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Technical Data

Terminal unit type	1L
Operating pressure for ejector drive gas	5 bar \pm 0.5 bar ¹⁾
Consumption – adjustable	10 to 30 L/min
Output	50 L/min
Temperature	+15 °C to +40 °C (+59 °F to +104 °F)
Classification	class I (vacuum)
As per EC Directive 93/42/EEC Annex IX	
UMDNS code	17-682
Universal Medical Device Nomenclature System (Nomenclature for medical devices)	

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1) 1 bar = 1 kPa x 100

2) Variable plaster thickness for wall mounting.

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
Directive 93/42/EEC
concerning Medical Devices

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90 52 657 - GA 6913.106 en/cs
© Drägerwerk AG & Co. KGaA
Edition: 1 – 2016-06
Dräger reserves the right to make
modifications to the equipment
without prior notice.

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