

Aesculap Endoscopic Technology

- GB** Instructions for use/Technical description
- USA** GN160 bipolar high-frequency surgery device
- D** Gebrauchsanweisung/Technische Beschreibung
Bipolares Hochfrequenz-Chirurgiegerät GN160
- F** Mode d'emploi/Description technique
Appareil de chirurgie bipolaire à haute fréquence GN160
- E** Instrucciones de manejo/Descripción técnica
Aparato de cirugía de alta frecuencia bipolar GN160
- I** Istruzioni per l'uso/Descrizione tecnica
Elettrobisturi ad alta frequenza bipolare GN160
- P** Instruções de utilização/Descrição técnica
Aparelho de electrocirurgia de alta frequência bipolar GN160
- NL** Gebruiksaanwijzing/Technische beschrijving
Bipolaire hoogfrequente chirurgische apparatuur GN160
- S** Bruksanvisning/Teknisk beskrivning
Bipolär högfrequens-kirurgiapparat GN160
- RUS** Инструкция по применению/Техническое описание
Биполярный высокочастотный хирургический прибор GN160
- CZ** Návod k použití/Technický popis
Bipolární vysokofrekvenční chirurgický přístroj GN160
- PL** Instrukcja użytkowania/Opis techniczny
Dwubiegunowe urządzenie chirurgiczne o wysokiej częstotliwości GN160
- SK** Návod na používanie/Technický opis
Bipolárne vysokofrekvenčné chirurgické zariadenie GN160
- TR** Kullanım Kılavuzu/Teknik açıklama
Bipolar yüksek frekanslı cerrahi cihazı GN160

Aesculap®

GN160 bipolar high-frequency surgery device

Legend

- 1 Power OFF switch
- 2 Power ON switch
- 3 Indicator lamp power ON
- 4 Operating mode display (STANDARD or FORCED)
- 5 Display field (STANDARD or MALIS value)
- 6 Unit display (STANDARD or MALIS)
- 7 Indicator lamp HF-ON
- 8 Output socket HF-bipolar
- 9 Indicator ring (green, white or red)
- 10 Control knob (with pushbutton function)
- 11 Connection socket (foot switch)
- 12 Equipotential connector
- 13 Mains power input socket
- 14 Fuse holder (with 2 fuses)
- 15 Nameplate

Symbols on product and packages

	Symbols: Caution
	Follow the instructions for use
	Foot switch
	Type CF unit with defibrillation protection
	Equipotentialization connector
	Fuse
	Alternating current
	Non-ionizing radiation
	Marking of electric and electronic devices according to directive 2002/96/EC (WEEE). see Disposal
	Manufacturer's batch designation
	Manufacturer's serial number
	Manufacturer's article number
	Date of manufacture

Contents

1.	Safe handling	
1.1	Safe for use in conformance with IEC regulations	
	The operating environment	
	Patient safety	
2.	Product description	
2.1	Scope of supply	
2.2	Components required for operation	
2.3	Intended use	
2.4	Operating principle	
2.5	Monitoring functions	
2.6	Output graphs	
2.7	Adaptation curve	
2.8	Maximum peak output voltage (UP)	
3.	Preparation and setup	
3.1	First use	
3.2	Stacking of units	
3.3	Moving stacks of units	
3.4	Dismantling stacks of devices	
4.	Working with the HF device GN160	
4.1	System set-up	
	Connecting the accessories	
	Connect power connection and switch on HF device GN160	
	Connecting a single-pedal foot switch (e.g. GN161)	
	Connecting bipolar accessories	
4.2	Function checks	
4.3	Safe operation	
	Switching the HF device on/off	
	Setting a STANDARD or MALIS value	
	Navigation in the submenu	
	Activating the HF current	
5.	Validated reprocessing procedure	
5.1	General safety instructions	
5.2	General information	
5.3	Preparation before cleaning	
5.4	Cleaning/disinfection	
	Product-specific safety instructions for the reprocessing procedure	
	Validated cleaning and disinfection procedure	
5.5	Wipe disinfection for electrical devices without sterilization	
5.6	Inspection, maintenance and checks	
5.7	Storage and transport	
6.	Maintenance	
6.1	Safety inspection	
7.	Troubleshooting list	
7.1	Fuse replacement	
8.	Technical Service	

9.	Accessories/Spare parts	20
10.	Technical data	20
10.1	Ambient conditions	21
11.	Disposal	21
12.	Distributor in the US/Contact in Canada for product information and complaints	21

1. Safe handling

CAUTION

Federal law restricts this device to sale by or on order of a physician!

Note

These instructions for use only describe the setup, function and operation of the GN160 high-frequency surgical device and is not suitable for initiating a newcomer to high-frequency surgery. A general description of the application of high-frequency electrosurgery can be found in the respective O.R. manuals and other relevant literature.

DANGER

Risk of injury to patients due to inappropriate application!

► **The product and accessories should only be operated by qualified or trained and experienced personnel.**

- Transport the product only in its original box.
- Prior to use, check that the product is in good working order.
- Observe "Notes on Electromagnetic Compatibility (EMC)", see TA022130.
- To prevent damage caused by improper setup or operation, and in order not to compromise warranty and manufacturers liability:
 - Use the product only in accordance with these instructions for use.
 - Follow the safety and maintenance instructions.
 - Only combine Aesculap products with each other.
 - Adhere to application instructions according to relevant norms, see Safe for use in conformance with IEC regulations.
- Inspect the accessories regularly: Electrode cables and endoscopic accessories, in particular, must be checked for possible damage to the insulation.
- Keep the instructions for use accessible for the user.
- Always adhere to applicable standards.

Aesculap®

GN160 bipolar high-frequency surgery device

1.1 Safe for use in conformance with IEC regulations

1.1.2 Patient safety

1.1.1 The operating environment

DANGER

Risk of death by electric shock!

- ▶ Do not open the product.
- ▶ Only ever connect the product to power mains with equipment grounding conductor.

WARNING

Danger of injury from burns and/or explosion from flammable gases! Sparks can occur even if the HF surgical device is used properly and according to its intended use.

- ▶ Do not use the device in explosion-hazard zones.
- ▶ When operating in the head or thoracic region, avoid using combustible anesthetics and accelerating gases (e.g. nitrous oxide or oxygen) or, when using such substances, ensure they are extracted from the region of operation.
- ▶ If possible, use incombustible cleaning and disinfecting agents.
- ▶ If combustible cleaning and disinfecting agents or solvents have to be used: Verify that such agents have evaporated prior to commencing HF surgery.
- ▶ Ensure that no combustible fluids accumulate under the patient or in body cavities (e.g. the vagina). Wipe off all fluids before using the HF surgical device.
- ▶ Ensure the absence of any endogenous, combustible gases.
- ▶ Check that oxygen-soaked materials (e.g. absorbent cotton or mull) are kept at a safe distance from the HF field, so that they cannot ignite.

CAUTION

Risk of interference with other devices!

Even during normal use, the HF surgical device creates electromagnetic fields that can interfere with other devices.

- ▶ Check that no electronic devices that could be subject to interference by electromagnetic fields are set up in the vicinity of the HF surgical device.

CAUTION

Restriction to view and/or side-effects due to the development of steam/smoke during HF surgery!

- ▶ If necessary, use smoke suction.

- ▶ Ensure that the device does not come in direct contact with the patient or in the sterile area respectively.
- ▶ Ensure that the user does not come into direct, simultaneous contact with the patient and the HF device.

DANGER

Danger due to inadequate preparation or faults in the HF surgical device!

- ▶ Check that the HF surgical device is in perfect working order.
- ▶ Ensure that neither the foot switch nor the hand switch has been penetrated by conductive fluids (e.g. blood, amniotic fluid).
- ▶ Ensure there is no short circuit in the foot or hand switch cables.

DANGER

Risk of burns suffered by the patient due to inadvertent activation of the HF surgical device!

- ▶ In case of any inadvertent activation of the HF surgical device, switch off the device immediately at the OFF/ON switch.
- ▶ Always exercise particular care when operating the foot switch.

DANGER

Risk of injury to the patient due to uncontrolled rise of the HF output voltage due to some fault in the HF surgical device!

- ▶ Stop using the HF surgical device as soon as it shows even the slightest anomaly.

WARNING

Risk of injury to patients/users due to defective power cord or missing protective ground connections!

- ▶ Check the mains power cord/protective ground connections.

WARNING

Danger of injuries due to muscle contraction, caused by stimulation of the nerves and muscles!

- ▶ Work with particular care on sensitive structures.

Always be sure to do the following:

- ▶ Position the patient in such a way that s/he is not in contact with any metal parts that are grounded or have a significant electric capacity against ground (e.g. operating table, fixtures). If necessary, interpose antistatic drape.
- ▶ Ensure that the patient will not be in contact with any damp cloths, drapes or bedding.
- ▶ Lay anti-static drape against parts of the body that tend to develop large amounts of sweat in order to safeguard against skin-to-skin contact with the torso.
- ▶ Siphon off urine with a catheter.
- ▶ For patients with cardiac pacemaker or other active implants, consult with the relevant medical specialist prior to applying HF surgery, so that irreparable damage to the pacemaker or implant can be avoided.
- ▶ If possible, remove from the vicinity of HF electrodes any electrodes from physiological monitoring devices that are not equipped with protective resistors or HF dampers.
- ▶ Do not use needle electrodes for intraoperative monitoring.
- ▶ Arrange the wires and cables of monitoring devices in such a way that they do not come in contact with the patient's skin.
- ▶ Keep the pathways of the wires and cables of HF electrodes as short as possible, and arrange them in such a way that they do not touch the patient or each other.
- ▶ Do not remove hot electrodes from the patient's body immediately after performing cutting or coagulation procedures.
- ▶ Never put down HF instruments or active electrodes on or next to the patient.
- ▶ Put down active electrodes that are not needed at any particular moment, in such a way that they will not touch the patient.
- ▶ Adjust the HF power output according to the intended surgical intervention, taking into account clinical experience and reference parameters.
- ▶ Always set the power output of HF devices to as low a level as possible.
- ▶ Should the output power appear insufficient with the usual settings, check that:
 - the working electrodes are clean,
 - all plug connections are properly in place.
- ▶ For operations involving unavoidable, continuous contact between the electrodes and the patient (e.g. endoscopic procedures), press the power OFF switch at the HF device immediately after any inadvertent activation of the electrode.
- ▶ Adjust the acoustic warning, which signals the activation of the electrode, to such a level that it will always be heard without difficulty.

2. Product description

2.1 Scope of supply

Designation	Art. no.
Bipolar HF surgical device	GN160
Notes on electromagnetic compatibility (EMC)	TA022130
GN160 bipolar high-frequency surgery device instructions for use	TA022247

2.2 Components required for operation

- Power cord
- Foot switch (single pedal)
- Bipolar connection cable
- Bipolar forceps

2.3 Intended use

The GN160 is a bipolar HF surgical device for coagulation with bipolar instruments (e.g. forceps). It can be used in all surgical specialties, also meeting the special requirements of neurosurgery (e.g. avoiding tissue adhesion to instruments, fine adjustment of HF power).

Aesculap®

GN160 bipolar high-frequency surgery device

2.4 Operating principle

The HF device GN160 is fitted with a bipolar output socket that can provide HF currents for coagulation (optimized for forceps).

The HF surgical device can be activated only through a single-pedal foot switch (e.g. footswitch for operation from any angle, GN161), which is connected via a socket in the rear panel of the HF surgical device.

The output power is set directly by turning control knob **10**, and is displayed in one of the two established units (STANDARD/MALIS; factory setting = STANDARD).

All other device parameters (e.g. operating mode, volume, etc.) and system information (software version, serial no, ...) are set and displayed via a submenu, which can be called up by pressing control knob **10**.

The operating mode display **4** and the display unit **6** for special functions are only visible when the corresponding special function is active. During normal operation, the display elements are hidden, contributing to the easy-to-use, uncluttered quality of the device controls.

The HF device GN160 offers two operating modes, both with a setting range of 1–60 STANDARD (or 10–170 MALIS):

- The STANDARD operating mode (factory default) allows rapid coagulation throughout the entire power range, without tissue adhesion or carbonization.
- The FORCED mode allows coagulation with increased sloughing of the surrounding tissue.

The HF device GN160 can also be combined with accessory devices (e.g. JET irrigation unit GN090).

2.5 Monitoring functions

Self-test

As soon as it is switched on, the HF unit performs a self-test of the control elements, the acoustic warning signal, the microprocessor and the hardware functions.

Continuous test cycle during operation

During operation, safety-relevant functions and signals are monitored through a continuous test cycle. As soon as a critical error is detected, the HF device terminates the HF activation.

The code number of the error is displayed in display field **5** and an acoustic warning signal may be issued, see Troubleshooting list.

2.6 Output graphs

Output power (P_{out}) as function of setting (STANDARD/MALIS):


Fig. 1 STANDARD


Fig. 2 MALIS

2.7 Adaptation curve

Operating modes STANDARD and FORCED: $(P_{out}) = f(R_L)$
(measured with original cable: $C_L = 175 \text{ pF}$)


Fig. 3 Operating mode: STANDARD


Fig. 4 Operating mode: FORCED

2.8 Maximum peak output voltage (U_p)

Note

The following graphs allow the user to judge whether the HF device or its output setting is suitable for a given accessory (insulation rating).

WARNING

Risk of injury to patients or users caused by inadequate accessories (insulation rating)!

► Make certain that the accessory voltage rating specified in the product documentation is higher than maximum peak output voltage for the intended setting.

Maximum peak output voltage (U_p) as function of setting (STANDARD/MALIS):


Fig. 5 STANDARD


Fig. 6 MALIS

Aesculap®

GN160 bipolar high-frequency surgery device

3. Preparation and setup

Non-compliance with the following instructions will preclude all responsibility and liability in this respect on the part of Aesculap.

- ▶ When setting up and operating the product, adhere to
 - national regulations for installation and operation,
 - national regulations on fire and explosion protection.

Note

For the safety of patients and users it is essential that the mains power cord and, especially, the protective earth connection are intact. In many cases defective or missing protective earth connections are not registered immediately.

- ▶ Connect the device via the potential equalization terminal at the rear panel of the device to the potential equalization system of the room used for medical purposes.

Note

The potential equalization lead can be ordered from the manufacturer as art. no. GK535 (4 m length) or TA008205 (0.8 m length).

3.1 First use

WARNING

Risk of injury and/or product malfunction due to incorrect operation of the electromedical system!

- ▶ Adhere to the instructions for use of any medical device.

3.2 Stacking of units

The stacking provisions integrated in the device cover allows secure mounting of an accessory device (e.g. JET Irrigation unit GN090) on top of the HF device GN160.

- ▶ Make certain the system is set up on a sufficiently stable support (e.g. a table, ceiling support, equipment cart, etc.).

Stacking products fitted with stacking devices:

- ▶ Do not exceed a stack height of 475 mm.
- ▶ Remove the caps from stacking cones **C**.
- ▶ Unscrew the feet of the product that will be on top.
- ▶ Move stacking cones **C** to the correct mounting position by applying a screwdriver onto the locking screw **B** and turning it counterclockwise to the positive stop.
- ▶ Place the upper device on the one below it.
- ▶ Securing the device: Turn locking screw **B** clockwise to the positive stop.
- ▶ Slightly lift the units to check that they are safely connected to each other.
- ▶ Place the units in a stable position.


Fig. 7 Stacking of units

Legend

- A Accessory device, e.g. JET irrigation device GN090
- B Locking screw
- C Stacking cone

3.3 Moving stacks of units

- ▶ Always lift stacks from underneath the bottommost unit.

3.4 Dismantling stacks of devices

Dismounting a stack of products fitted with stacking devices:

- ▶ Apply a screwdriver on locking screw **B** and turn it counterclockwise to the positive stop.
- ▶ Remove the topmost device.
- ▶ Screw back on the feet of the device, either by hand or using a screwdriver.


Fig. 8 Dismantling stacks of devices

Legend

B Locking screw

4. Working with the HF device GN160

4.1 System set-up

4.1.1 Connecting the accessories

Risk of injury due to unapproved configuration using additional components!

- ▶ For all applied components, ensure that their classification matches that of the application component (e.g. Type BF or Type CF) of the respective device.

Combinations of accessories that are not mentioned in the present instructions for use may only be employed if they are specifically intended for the respective application, and if they do not compromise the performance and safety characteristics of the products.

Also note that any equipment connected at the interfaces must demonstrably meet the respective IEC standards (e.g. IEC 60950 for data processing equipment, IEC/DIN EN 60601-1 for electromedical devices).

All configurations must comply with basic standard IEC/DIN EN 60601-1. Any individual connecting devices with one another is responsible for such configuration and must ensure compliance with basic standard IEC/DIN EN 60601-1 or applicable national standards.

The insulation of the accessories (e.g. HF cables, instruments) must be adequately rated for the maximum peak output voltage (see IEC DIN EN 60601-2-2).

- ▶ Please address your B. Braun/Aesculap partner or Aesculap Technical Service with any inquiries in this respect; for a contact address, see Technical Service.

4.1.2 Connect power connection and switch on HF device GN160

Risk of fatal injury from electric shock!

- ▶ Connect the product only to a grounded power supply.

The mains voltage must correspond to the voltage indicated on the type plate at the back of the unit.

The HF device is fitted with a universal power adapter for mains voltages 100–120 V and 220–240 V.

- ▶ Plug in the power cord in mains power input socket **13** in the rear panel of the HF device.
- ▶ Connect power supply by inserting into a wall socket.
- ▶ Switch on the HF device with the "Power ON" switch **2**.

Indicator lamp power ON **3** is illuminated.

The HF device performs a self-test.

On switching on the device for the first time, the basic factory default settings will be shown on the display at the end of the self-test, see Fig. 9.


Fig. 9 Factory default setting

On every further switch-on, the device settings last used prior to the most recent switch-off will be applied and displayed.

Aesculap®

GN160 bipolar high-frequency surgery device

4.1.3 Connecting a single-pedal foot switch (e.g. GN161)

Both the foot switch plug and connection socket (foot switch) 11 carry an arrow marking. To connect the plug to the connection socket correctly, the two arrows must be aligned.

Fig. 10 Connecting a single-pedal foot switch (e.g. GN161)

- Plug in the foot switch plug in connection socket foot switch 11.

4.1.4 Connecting bipolar accessories

WARNING

Risk for patients and users due to inadvertent current supply to an instrument!

- **Make certain that not more than one patient cable is connected to the output socket.**

Note

Connect the bipolar cable only after switching on the GN160!

The following plug types can be connected at output socket HF 8.


Fig. 11 Coax plug


Fig. 12 2-pin plug

Note

IEC/DIN EN 60601-2-2:2007 requires active plugs with more than one contact pin to be specified with fixed pin spacing. Loose (flying) leads are not allowed!

- Plug in the bipolar cable in output socket HF 8.

4.2 Function checks

Factory settings display:

Operating mode: Standard

Unit: STANDARD

Note

Every subsequent switch-on of the HF device is followed by the display of the last device settings applied prior to the most recent switch-off. Any changed settings applied during the function test are saved immediately.

Any device settings other than the factory settings must be noted down before the function test is continued!

If in doubt, we recommend restoring the factory settings.

- ▶ Always carry out function checks before using the product.
- ▶ Check that the accessories do not show any visible damage.
- ▶ Prepare and set up the HF unit, see Preparation and setup.
- ▶ Connect the unit to power supply, see Connect power connection and switch on HF device GN160.
- ▶ Check the functionality of the following elements, one after the other, in the following sequence:

User action	Device reaction
1. Press power ON switch 2 .	Indicator lamp power ON 3 is illuminated. The HF device performs a self-test, accompanied by an acoustic signal. All display elements are activated in rapid sequence. Sequence: 1. Display operating mode (FORCED) 4 2. All segments of display field 5 (as running light) 3. Display Unit (MALIS) 6 4. Indicator lamp HF-ON 7 5. Indicator ring (green) 9 6. Indicator ring (white) 9 7. Indicator ring (red) 9 An acoustic signal confirms that the self-test is complete and the device ready for operation. Indicator ring 9 is illuminated green after completion of the self-test. The most recent device settings applied prior to the last switch-off are displayed.
2. Rotate the adjustment knob 10 in a clockwise direction.	The parameter value displayed in display field 5 is increased.

User action	Device reaction
3. Rotate the adjustment knob 10 in a counterclockwise direction.	The parameter value displayed in display field 5 is decreased.
4. Press control knob 10 (>2 s).	Indicator ring 9 is illuminated white.  Display of selected operating mode in display field 5 Operating mode STANDARD  Operating mode FORCED 
5. Press control knob 10 again (>2 s).	The submenu is quitted. Indicator ring 9 is illuminated green.  Display of most recent setting in display field 5
6. Plug in the single-pedal foot switch plug in connection socket (foot switch) 11 .	No device reaction If the HF device is activated at this stage, the foot switch is defective!
7. Press the single-pedal foot switch.	The HF device is activated. Indicator ring 9 turns off. Indicator lamp HF-ON 7 illuminated If there is no device reaction at this stage, there can be a fault either in the HF device or in the foot switch. Repeat the test with another foot switch, see Troubleshooting list.
8. Press power OFF switch 1 .	Indicator lamp power ON 3 turns off.

- ▶ In case of device reactions other than those intended: Do not use the HF device and contact Technical Service, see Technical Service.

Aesculap®

GN160 bipolar high-frequency surgery device

4.3 Safe operation

Risk of injury and/or malfunction!
 ► Always carry out a function check prior to using the product.

Risk of injury when using the product beyond the field of view!
 ► Apply the product only under visual control.

4.3.1 Switching the HF device on/off

- Switching on the HF device: Press power ON switch 2.
- Switching off the HF device: Press power OFF switch 1.

4.3.2 Setting a STANDARD or MALIS value

The STANDARD or MALIS value is set via control knob 10 and displayed in display field 5.

Slow rotation of control knob 10 allows selecting any possible setting; quick rotation allows going through the entire setting range by only a few turns.

For both operating types (STANDARD/FORCED), STANDARD values of 1-60 and MALIS values of 10-170 respectively are available.

The possible setting range and graduation depend on the unit selected (STANDARD or MALIS):

Display unit	Setting range	Graduation
STANDARD	1,0–10,0	0,1
	10,0–30,0	0,5
	30–60	1
MALIS	10–170	1

Note

The precision of the setting meets the requirements of IEC/DIN EN 60601-2-2!

Note

At the start of the setting process, the HF unit issues a short acoustic signal, alerting the user to the impending change of the STANDARD or MALIS setting. In this way it is ensured that the user will be aware of any inadvertent change of the setting!

- Turn control knob 10 clockwise or counterclockwise to set the intended STANDARD or MALIS value.

4.3.3 Navigation in the submenu

The submenu allows the following settings and display of system information:

- Settings
 - Operating mode (STANDARD/FORCED)
 - Noise
 - Unit (STANDARD/MALIS)
- System information
 - Last error event
 - Software version
 - Serial number

Note

The HF device can not be activated while the submenu is active!

- To call up the submenu: Press and hold control knob 10 (>2 s). The first submenu item (Operating mode) is displayed.
- To call up the next submenu item: Briefly press control knob 10 (<2 s). The next submenu item is called up.
- To quit the submenu: Press and hold control knob 10 (>2 s). The setting for the displayed submenu item is applied. If no input via the control knob is applied through a period of 30 s, the setting of the displayed submenu item is applied and the submenu is quitted automatically.

Setting the operating mode

- ▶ To select an operating mode: Press and hold control knob **10** (>2 s).

The following display appears in display field 5:

- In operating mode STANDARD

- In operating mode FORCED

- ▶ To set an operating mode: Turn control knob **10**.

The following display appears in display field 5:

- In operating mode STANDARD

Display Operating mode **4** is not visible.

The following display appears in display field 5:

- In operating mode FORCED

Display Operating mode **4** is illuminated.

- ▶ To apply the setting: Briefly press control knob **10** (<2 s).
The next submenu item is called up.
- ▶ To quit the submenu: Press and hold control knob **10** (>2 s).
The setting for the displayed submenu item is applied.

Adjusting the volume

The signal volume of the HF device can be adjusted to any of 9 volume levels.

- ▶ To select an operating mode: Press and hold control knob **10** (>2 s).
The first submenu item (Operating mode) is displayed.
- ▶ Repeatedly, briefly (<2 s) press control knob **10** until submenu item (Volume) is displayed.

The following display appears in display field 5 (SP (for speaker) and the actual volume level setting (1-9)):

An acoustic signal is issued at the currently set volume.

- ▶ To adjust the volume: Turn control knob **10**.
The actual volume level is displayed and an acoustic signal is issued at the actual volume.
- ▶ To adopt the setting: Briefly press control knob **10** (<2 s).
The next submenu item is called up.

- ▶ To quit the submenu: Briefly press control knob **10** (<2 s).
The setting for the displayed submenu item is applied.
- or -
- ▶ To quit the submenu without adopting the setting: Press control knob **10** for some time (>2 s) or wait for 30 s.

Selecting the unit (STANDARD/MALIS)

- ▶ To select the unit: Press and hold control knob **10** (>2 s).
The submenu item (Operating mode) is displayed.
- ▶ Repeatedly, briefly (<2 s) press control knob **10** until submenu item Unit (Unt) (Unt = Unit) is displayed.

As long as the unit STANDARD is selected, the Unit display **6** is hidden. If the unit MALIS is selected, display Unit **6** is illuminated.

- ▶ To select the unit: Turn control knob **10**.
- ▶ To adopt the setting: Briefly press control knob **10** (<2 s).
The next submenu item is called up.
- ▶ To quit the submenu: Press control knob **10** briefly (<2 s) or wait for 30 s.
The setting for the displayed submenu item is applied.

Displaying the last error event

Note

For more detailed fault descriptions, see *Troubleshooting list*.

- ▶ To call up the submenu: Press and hold control knob **10** (>2 s).
The first submenu item (Operating mode) is displayed.
- ▶ Repeatedly, briefly (<2 s) press control knob **10** until submenu item (Last error event) is displayed.
The following display appears in display field 5 (E (for error) and the respective error code):

- ▶ To switch to another submenu item: Briefly press control knob **10** (<2 s).
The next submenu item is called up.
- ▶ To quit the submenu: Press control knob **10** briefly (<2 s) or wait for 30 s.

Aesculap®

GN160 bipolar high-frequency surgery device

Displaying the software version

- ▶ Press and hold control knob **10** (>2 s).
The submenu item (Operating mode) is displayed.
- ▶ Repeatedly, briefly (<2 s) press control knob **10** until submenu item (Software version) is displayed.
The following marquee display appears in display field **5** (SoF (for software), followed by an underscore and the three-digit version number):

- ▶ To switch to another submenu item: Briefly press control knob **10** (<2 s).
The next submenu item is called up.
- ▶ To quit the submenu: Press control knob **10** briefly (<2 s) or wait for 30 s.

Displaying the serial number

- ▶ Press and hold control knob **10** (>2 s).
The submenu item (Operating mode) is displayed.
- ▶ Repeatedly, briefly (<2 s) press control knob **10** until submenu item (Serial number) is displayed.
The following marquee display appears in display field **5** (Sn for serial number, followed by an underscore and the four-digit device serial number):

- ▶ To quit the submenu: Briefly press control knob **10** (<2 s).
- or -
- ▶ Press control knob **10** for some time (>2 s) or wait for 30 s.

4.3.4 Activating the HF current

- ▶ Be sure that the patient is prepared in such a way that the HF unit can be applied without any risk to the patient.
- ▶ Verify that the required accessories are properly connected and plugged in.
- ▶ Be sure that all settings on the HF device are appropriate for the intended application.
- ▶ Press the foot switch.
Indicator lamp HF-ON **7** will be illuminated for the duration of the activation.

Note

The activation period is limited to a maximum of 120 s, meaning: after 120 s of continuous activation, the HF current is deactivated automatically. To alert the user to the automatic deactivation, the activation signal volume will rise, during the last 60 s, from the selected volume setting to maximum volume.

5. Validated reprocessing procedure

5.1 General safety instructions

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

The recommended chemistry was used for validation.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For the latest information on reprocessing and material compatibility see also the Aesculap extranet at www.extranet.bbraun.com

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

5.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

5.4.2 Validated cleaning and disinfection procedure

Validated procedure	Special features	Reference
Wipe disinfection for electrical devices without sterilization	–	Chapter Wipe disinfection for electrical devices without sterilization

5.3 Preparation before cleaning

- ▶ Keep the product separate from the power supply.
- ▶ Remove accessories.

5.4 Cleaning/disinfection

5.4.1 Product-specific safety instructions for the reprocessing procedure

Risk of electric shock and fire hazard!

- ▶ Unplug the device before cleaning.
- ▶ Do not use flammable or explosive cleaning or disinfecting solutions.
- ▶ Ensure that no fluids will penetrate the product.

Damage to, or destruction of the product caused by mechanical cleaning/disinfection!

- ▶ Only clean and disinfect the product manually.
- ▶ Do not sterilize the product under any circumstances.

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Only use cleaning/disinfecting agents approved for surface cleaning. Follow the manufacturer's instructions for the respective cleaning/disinfecting agent.
- ▶ Observe specifications regarding concentration, temperature and exposure time.

Aesculap®

GN160 bipolar high-frequency surgery device

5.5 Wipe disinfection for electrical devices without sterilization

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Wipe disinfection	RT	≥1	-	-	Meliseptol HBV wipes 50 % Propan-1-ol

RT: Room temperature

Phase I

- ▶ Remove any visible residues with a disposable disinfectant wipe.
- ▶ Wipe all surfaces of the optically clean product with a fresh, disposable disinfectant wipe.
- ▶ Observe the specified application time (1 min minimum).

5.6 Inspection, maintenance and checks

- ▶ Inspect the product after each cleaning and disinfecting cycle to be sure it is: clean, functional, and undamaged.
- ▶ Inspect all cables, and in particular electrode cables, for any damage to their insulation.
- ▶ Set aside the product if it is damaged.
- ▶ Perform a safety inspection annually, see Maintenance.

5.7 Storage and transport

Damage to the product caused by switch-on to early after storage/transport at temperature below +10 °C (50° F)!

- ▶ **Allow the HF device to acclimate at room temperature for about 3 hours.**

- ▶ Transport the product only in its original box.
- ▶ Observe storage and transport conditions, see Ambient conditions.

6. Maintenance

No servicing is required for the HF device GN160.

Note

Maintenance must only be carried out by authorized personnel. If necessary, the circuit diagrams and the service manual, which contains all the necessary documents, are provided for this.

6.1 Safety inspection

A safety inspection must be carried out annually.

The person in charge of the inspection will document the inspection findings and measurement values according to the printed inspection protocol.

- ▶ Only have the product and its accessories inspected by persons who possess the requisite training, knowledge and experience and who are authorized to do so.
- ▶ If there are significant deviations from the inspection protocol values, or if the prescribed maximum values are exceeded: Send in device.

For technical service, please contact your national B. Braun/Aesculap agency, see Technical Service.

		Inspection protocol – safety inspection				
TEST INTERVAL: 1 Year TYPE OF DEVICE: Bipolar Electrosurgical Unit GN160 MANUFACTURED BY: Aesculap AG, Am Aesculap-Platz 78532 Tuttlingen/Germany		RESPONSIBLE SN: INVENTORY NO.: ORGANIZATION:				
SCOPE OF INSPECTION					Finding	
1.) VISUAL INSPECTION					PASS	FAIL
1.1 Inspection of power cable for visually detectable damage						
1.2 Inspection of foot control for visually detectable damage						
1.3 Inspection of coagulation cable for visible signs of damage (record their number) and check the fitting of the plug connector						
1.4 Type plate still present and readable						
1.5 Check fuse links for prescribed values						
1.6 Overall condition of device: dirt, damage						
2.) ELEKTRICAL INSPECTION					Measured value	OK Yes / No
2.1 Protective earth resistance incl. mains cable according IEC62353:2007				$R_{max} = 0,3 \Omega$ (at $\geq 0,2 A$)		
2.2 Earth leakage current according IEC60601-1*		N.C. $I_{max} = 0,50 mA$ S.F.C. $I_{max} = 1,00 mA$				
2.3 Enclosure leakage current/touch current according IEC60601-1*		N.C. $I_{max} = 0,10 mA$ S.F.C. $I_{max} = 0,50 mA$				
2.4 Patient leakage current according IEC60601-1*		N.C. $I_{max} = 0,01 mA$ S.F.C. $I_{max} = 0,05 mA$				
2.5 Equipment leakage current - alternative method according IEC62353:2007				$I_{max} = 1,00 mA$		
2.6 Applied part leakage current - alternative method according IEC62353:2007				$I_{max} = 0,05 mA$		
2.7 HF-leakage current bipolar according IEC60601-2-2:2009				$I_{max} = 55 mA$		
2.8 Insulation resistance power according IEC62353:2007 (Preparation according Service Manual)				test voltage = 500 V DC		
2.8.1 Mains against HF output				$R_{min} = 7 M\Omega$		
2.8.2 Mains against housing				$R_{min} = 2 M\Omega$		
2.8.3 Bipolar HF-output against housing				$R_{min} = 2 M\Omega$		
2.9 HF-output power measurement with induction-free resistance						
HF-output		Operational mode	Dosage	R_{Load}	Rated value	
Bipolar		STANDART	60	75 Ω	60 Watt $\pm 20 \%$	
		FORCED	60	100 Ω	60 Watt $\pm 20 \%$	
2.10 DC resistance between the two HF output poles according IEC60601-2-2:2009 section 201.8.4.102				$R_{min} = 2 M\Omega$		
2.11 Function test according to instructions for use					Performed:	
Place/Date		Inspector/Signatur		Responsible Organization		

* The listed allowable values comply with the requirements of the IEC60601-1:1988+A1:1991+A2:1995 and IEC60601-1:2005

Fig. 13 Test protocol

Aesculap®

GN160 bipolar high-frequency surgery device

7. Troubleshooting list

If an error occurs, the following display appears in display field 5, which is made up of the letters E (for error) and the error code:

At the same time, the indicator ring 9 lights up red.

Mal-function	Detection	Cause	Remedy
1	Switch-on self-test activated for foot switch	Foot switch pressed during switch-on, or foot switch defective	Release or unplug the foot switch Aesculap Contact Technical Service, see Technical Service
2	Switch-on self-test control knob	Control knob pressed during switch-on, or control knob defective	Release the control knob Aesculap Contact Technical Service, see Technical Service
3	Control knob pushbutton function	Control knob kept pressed for too long; control knob defective	Release the control knob Aesculap Contact Technical Service, see Technical Service
4	Time limit exceeded	HF activation time limit exceeded	Release the foot switch
10-89	Internal error	Internal device fault/error	Switch off the HF unit and switch it on again Aesculap Contact Technical Service, see Technical Service
90-99	Calibration	Lack of balancing	Aesculap Contact Technical Service, see Technical Service

7.1 Fuse replacement

Risk of fatal injury from electric shock!
 ► Unplug the device before changing the fuses!

Specified fuse sets	TA020005(2 pcs.)
Designation	T 2.5 AH / 250 V
Time-Current characteristic	T (idle)
Breaking capacity	H (1 500 A)
Construction	5 mm x 20 mm

- Push down latch nose **A** at fuse holder **C** and unlock the fuse holder.
- Remove fuse holder **C**.
- Replace both fuse sets **B**.
- Insert fuse holder **C** so that it clicks into its proper position.


Fig. 14 Fuse replacement

Legend

- A** Latch nose
- B** Fuse sets
- C** Fuse holder

Note

If the fuses burn out frequently, the device is faulty and should be repaired, see Technical Service.

8. Technical Service

Risk of injury and/or malfunction!
 ► Do not modify the product.

- For service and repairs, please contact your national B. Braun/ Aesculapagency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

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Phone: +1 (800) 214 -3392

Fax: +1 (314) 895 -4420

Other service addresses can be obtained from the address indicated above.

Aesculap®

GN160 bipolar high-frequency surgery device

9. Accessories/Spare parts

Power cord

Art. no.	Approvals	Length [m]
TE730	Europe	5.0
TE734	United Kingdom and Ireland	5.0
TE735	USA, Canada, Japan	3.5

Potential-compensation line

Art. no.	Designation	Length [m]
GK535	Potential-compensation line	4.0
TA008205	Potential-compensation line	0.8

Foot switch

Art. no.	Designation
GN161	Single-pedal foot switch (round)
GK226	Single foot switch

Bipolar instruments and connection cable

Note

Information concerning the instruments and connection cable is available on request and can be found in the Aesculap brochure C-304-81.

Spare parts

Art. no.	Designation
TA020005	Fuses that can be used with this unit
TA022130	Notes on electromagnetic compatibility (EMC)

Note

Information on other accessories and replacement parts is available on request and can be found in Aesculap brochure C-304-81.







10. Technical data

Classification acc. to Directive 93/42/EEC

Art. no.	Designation	Class
GN160	Bipolar HF surgical device	IIb

Mains voltage ranges	100-120 V ± 10 % 220-240 V ± 10 %
Frequency	50-60 Hz
Power consumption (with a 60-Watt HF delivery)	1.00-0.81 A (for 100-120 V) 0.44-0.41 A (for 220-240 V)
Power consumption (service-ability)	ca. 0,12 A (for 100-120 V) ca. 0.14 A (for 220-240 V)
Protection class (acc. to IEC/ DIN EN 60601-1)	I
Protection type acc. to IEC 60529	IP31
Foot switch circuit	fused in conformance with IEC 60601, certified for use in a "medical environment"
Output power	STANDARD: 60 W at 75 Ohms FORCED: 60 W at 100 Ohms
Output frequency	346 kHz
Operating mode	Int 10 s/30 s
Hardware protection	T 2.5 AH / 250 V
Time-Current characteristic	T (idle)
Breaking capacity	H (1 500 A)
Construction	5 mm x 20 mm
Weight	6 kg
Dimensions (L x W x H)	325 mm x 305 mm x 135 mm
Applied part	Type CF Defibrillation-protected output Recovery time 0 s
EMC	IEC/DIN EN 60601-1-2 CISPR11 Class A
Conforming to standard	IEC/DIN EN 60601-1 IEC/DIN EN 60601-2-2

10.1 Ambient conditions

	Operation	Storage and transport
Temperature		
Relative humidity		
Atmospheric pressure		

Note

The atmospheric pressure of 700 hPa corresponds to a maximum operating altitude of 3 000 m.

11. Disposal

Note

The user institution is obliged to process the product before its disposal, see Validated reprocessing procedure.

 	<p>Adhere to national regulations when disposing of or recycling the product, its components and its packaging! The recycling pass can be downloaded from the Extranet as a PDF document under the respective article number. (The recycling pass includes disassembling instructions for the product, as well as information for proper disposal of components harmful to the environment.)</p> <p>Products carrying this symbol are subject to separate collection of electrical and electronic devices. Within the European Union, disposal is taken care of by the manufacturer as a free-of-charge service.</p>
--	--

- Detailed information concerning the disposal of the product is available through your national B. Braun/Aesculap agency, see Technical Service.

12. Distributor in the US/Contact in Canada for product information and complaints

3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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		Inspection protocol – safety inspection				
TEST INTERVAL: 1 Year TYPE OF DEVICE: Bipolar Electrosurgical Unit GN160 MANUFACTURED BY: Aesculap AG, Am Aesculap-Platz 78532 Tuttlingen/Germany		RESPONSIBLE SN: INVENTORY NO.: ORGANIZATION:				
SCOPE OF INSPECTION					Finding	
1.) VISUAL INSPECTION					PASS	FAIL
1.1	Inspection of power cable for visually detectable damage					
1.2	Inspection of foot control for visually detectable damage					
1.3	Inspection of coagulation cable for visible signs of damage (record their number) and check the fitting of the plug connector					
1.4	Type plate still present and readable					
1.5	Check fuse links for prescribed values					
1.6	Overall condition of device: dirt, damage					
2.) ELEKTRICAL INSPECTION					Measured value	OK Yes / No
2.1	Protective earth resistance incl. mains cable according IEC62353:2007		R _{max} = 0,3 Ω (at ≥0,2 A)			
2.2	Earth leakage current according IEC60601-1*		N.C. S.F.C.	I _{max} = 0,50 mA I _{max} = 1,00 mA		
2.3	Enclosure leakage current/touch current according IEC60601-1*		N.C. S.F.C.	I _{max} = 0,10 mA I _{max} = 0,50 mA		
2.4	Patient leakage current according IEC60601-1*		N.C. S.F.C.	I _{max} = 0,01 mA I _{max} = 0,05 mA		
2.5	Equipment leakage current - alternative method according IEC62353:2007		I _{max} = 1,00 mA			
2.6	Applied part leakage current - alternative method according IEC62353:2007		I _{max} = 0,05 mA			
2.7	HF-leakage current bipolar according IEC60601-2-2:2009		I _{max} = 55 mA			
2.8	Insulation resistance power according IEC62353:2007 (Preparation according Service Manual)		test voltage = 500 V DC			
	2.8.1	Mains against HF output		R _{min} = 7 MΩ		
	2.8.2	Mains against housing		R _{min} = 2 MΩ		
	2.8.3	Bipolar HF-output against housing		R _{min} = 2 MΩ		
2.9 HF-output power measurement with induction-free resistance						
	HF-output	Operational mode	Dosage	R _{Load}	Rated value	
	Bipolar	STANDART	60	75 Ω	60 Watt ± 20 %	
		FORCED	60	100 Ω	60 Watt ± 20 %	
2.10	DC resistance between the two HF output poles according IEC60601-2-2:2009 section 201.8.4.102		R _{min} = 2 MΩ			
2.11	Function test according to instructions for use				Performed:	
Place/Date		Inspector/Signatur		Responsible Organization		

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		Inspection protocol – safety inspection				
TEST INTERVAL: 1 Year TYPE OF DEVICE: Bipolar Electrosurgical Unit GN160 MANUFACTURED BY: Aesculap AG, Am Aesculap-Platz 78532 Tuttlingen/Germany		RESPONSIBLE SN: INVENTORY NO.: ORGANIZATION:				
SCOPE OF INSPECTION					Finding	
1.) VISUAL INSPECTION					PASS	FAIL
1.1	Inspection of power cable for visually detectable damage					
1.2	Inspection of foot control for visually detectable damage					
1.3	Inspection of coagulation cable for visible signs of damage (record their number) and check the fitting of the plug connector					
1.4	Type plate still present and readable					
1.5	Check fuse links for prescribed values					
1.6	Overall condition of device: dirt, damage					
2.) ELEKTRICAL INSPECTION					Measured value	OK Yes / No
2.1	Protective earth resistance incl. mains cable according IEC62353:2007		R _{max} = 0,3 Ω (at ≥0,2 A)			
2.2	Earth leakage current according IEC60601-1*		N.C. S.F.C.	I _{max} = 0,50 mA I _{max} = 1,00 mA		
2.3	Enclosure leakage current/touch current according IEC60601-1*		N.C. S.F.C.	I _{max} = 0,10 mA I _{max} = 0,50 mA		
2.4	Patient leakage current according IEC60601-1*		N.C. S.F.C.	I _{max} = 0,01 mA I _{max} = 0,05 mA		
2.5	Equipment leakage current - alternative method according IEC62353:2007		I _{max} = 1,00 mA			
2.6	Applied part leakage current - alternative method according IEC62353:2007		I _{max} = 0,05 mA			
2.7	HF-leakage current bipolar according IEC60601-2-2:2009		I _{max} = 55 mA			
2.8	Insulation resistance power according IEC62353:2007 (Preparation according Service Manual)		test voltage = 500 V DC			
	2.8.1	Mains against HF output	R _{min} = 7 MΩ			
	2.8.2	Mains against housing	R _{min} = 2 MΩ			
	2.8.3	Bipolar HF-output against housing	R _{min} = 2 MΩ			
2.9 HF-output power measurement with induction-free resistance						
	HF-output	Operational mode	Dosage	R _{Load}	Rated value	
	Bipolar	STANDART	60	75 Ω	60 Watt ± 20 %	
		FORCED	60	100 Ω	60 Watt ± 20 %	
2.10	DC resistance between the two HF output poles according IEC60601-2-2:2009 section 201.8.4.102		R _{min} = 2 MΩ			
2.11	Function test according to instructions for use				Performed:	
Place/Date		Inspector/Signatur		Responsible Organization		

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B BRAUN SHARING EXPERTISE		Inspection protocol – safety inspection				
TEST INTERVAL:		1 Year				
TYPE OF DEVICE:		Bipolar Electrosurgical Unit GN160				
MANUFACTURED BY:		Aesculap AG, Am Aesculap-Platz 78532 Tuttlingen/Germany				
SN:		INVENTORY NO.:		RESPONSIBLE ORGANIZATION:		
SCOPE OF INSPECTION					Finding	
1.) VISUAL INSPECTION					PASS	FAIL
1.1	Inspection of power cable for visually detectable damage					
1.2	Inspection of foot control for visually detectable damage					
1.3	Inspection of coagulation cable for visible signs of damage (record their number) and check the fitting of the plug connector					
1.4	Type plate still present and readable					
1.5	Check fuse links for prescribed values					
1.6	Overall condition of device: dirt, damage					
2.) ELEKTRICAL INSPECTION					Measured value	OK Yes / No
2.1	Protective earth resistance incl. mains cable according IEC62353:2007			$R_{max} = 0,3 \Omega$ (at $\geq 0,2 A$)		
2.2	Earth leakage current according IEC60601-1*			N.C. S.F.C.	$I_{max} = 0,50 mA$ $I_{max} = 1,00 mA$	
2.3	Enclosure leakage current/touch current according IEC60601-1*			N.C. S.F.C.	$I_{max} = 0,10 mA$ $I_{max} = 0,50 mA$	
2.4	Patient leakage current according IEC60601-1*			N.C. S.F.C.	$I_{max} = 0,01 mA$ $I_{max} = 0,05 mA$	
2.5	Equipment leakage current - alternative method according IEC62353:2007			$I_{max} = 1,00 mA$		
2.6	Applied part leakage current - alternative method according IEC62353:2007			$I_{max} = 0,05 mA$		
2.7	HF-leakage current bipolar according IEC60601-2-2:2009			$I_{max} = 55 mA$		
2.8	Insulation resistance power according IEC62353:2007 (Preparation according Service Manual)			test voltage = 500 V DC		
	2.8.1	Mains against HF output		$R_{min} = 7 M\Omega$		
	2.8.2	Mains against housing		$R_{min} = 2 M\Omega$		
	2.8.3	Bipolar HF-output against housing		$R_{min} = 2 M\Omega$		
	2.9 HF-output power measurement with induction-free resistance					
	HF-output	Operational mode	Dosage	R_{Load}	Rated value	
	Bipolar	STANDART	60	75 Ω	60 Watt $\pm 20 \%$	
		FORCED	60	100 Ω	60 Watt $\pm 20 \%$	
	2.10 DC resistance between the two HF output poles according IEC60601-2-2:2009 section 201.8.4.102			$R_{min} = 2 M\Omega$		
	2.11 Function test according to instructions for use				Performed:	
Place/Date		Inspector/Signatur		Responsible Organization		

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		Inspection protocol – safety inspection				
TEST INTERVAL: 1 Year TYPE OF DEVICE: Bipolar Electrosurgical Unit GN160 MANUFACTURED BY: Aesculap AG, Am Aesculap-Platz 78532 Tuttlingen/Germany		RESPONSIBLE SN: INVENTORY NO.: ORGANIZATION:				
SCOPE OF INSPECTION					Finding	
1.) VISUAL INSPECTION					PASS	FAIL
1.1 Inspection of power cable for visually detectable damage						
1.2 Inspection of foot control for visually detectable damage						
1.3 Inspection of coagulation cable for visible signs of damage (record their number) and check the fitting of the plug connector						
1.4 Type plate still present and readable						
1.5 Check fuse links for prescribed values						
1.6 Overall condition of device: dirt, damage						
2.) ELEKTRICAL INSPECTION					Measured value	OK Yes / No
2.1 Protective earth resistance incl. mains cable according IEC62353:2007				$R_{max} = 0,3 \Omega$ (at $\geq 0,2 A$)		
2.2 Earth leakage current according IEC60601-1*		N.C. $I_{max} = 0,50 mA$ S.F.C. $I_{max} = 1,00 mA$				
2.3 Enclosure leakage current/touch current according IEC60601-1*		N.C. $I_{max} = 0,10 mA$ S.F.C. $I_{max} = 0,50 mA$				
2.4 Patient leakage current according IEC60601-1*		N.C. $I_{max} = 0,01 mA$ S.F.C. $I_{max} = 0,05 mA$				
2.5 Equipment leakage current - alternative method according IEC62353:2007				$I_{max} = 1,00 mA$		
2.6 Applied part leakage current - alternative method according IEC62353:2007				$I_{max} = 0,05 mA$		
2.7 HF-leakage current bipolar according IEC60601-2-2:2009				$I_{max} = 55 mA$		
2.8 Insulation resistance power according IEC62353:2007 (Preparation according Service Manual)				test voltage = 500 V DC		
2.8.1 Mains against HF output				$R_{min} = 7 M\Omega$		
2.8.2 Mains against housing				$R_{min} = 2 M\Omega$		
2.8.3 Bipolar HF-output against housing				$R_{min} = 2 M\Omega$		
2.9 HF-output power measurement with induction-free resistance						
HF-output		Operational mode	Dosage	R_{load}	Rated value	
Bipolar		STANDART	60	75 Ω	60 Watt $\pm 20 \%$	
		FORCED	60	100 Ω	60 Watt $\pm 20 \%$	
2.10 DC resistance between the two HF output poles according IEC60601-2-2:2009 section 201.8.4.102				$R_{min} = 2 M\Omega$		
2.11 Function test according to instructions for use					Performed:	
Place/Date		Inspector/Signatur		Responsible Organization		

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		Inspection protocol – safety inspection				
TEST INTERVAL: 1 Year TYPE OF DEVICE: Bipolar Electrosurgical Unit GN160 MANUFACTURED BY: Aesculap AG, Am Aesculap-Platz 78532 Tuttlingen/Germany		RESPONSIBLE SN: INVENTORY NO.: ORGANIZATION:				
SCOPE OF INSPECTION					Finding	
1.) VISUAL INSPECTION					PASS	FAIL
1.1 Inspection of power cable for visually detectable damage						
1.2 Inspection of foot control for visually detectable damage						
1.3 Inspection of coagulation cable for visible signs of damage (record their number) and check the fitting of the plug connector						
1.4 Type plate still present and readable						
1.5 Check fuse links for prescribed values						
1.6 Overall condition of device: dirt, damage						
2.) ELEKTRICAL INSPECTION					Measured value	OK Yes / No
2.1 Protective earth resistance incl. mains cable according IEC62353:2007				$R_{max} = 0,3 \Omega$ (at $\geq 0,2$ A)		
2.2 Earth leakage current according IEC60601-1*		N.C. $I_{max} = 0,50$ mA S.F.C. $I_{max} = 1,00$ mA				
2.3 Enclosure leakage current/touch current according IEC60601-1*		N.C. $I_{max} = 0,10$ mA S.F.C. $I_{max} = 0,50$ mA				
2.4 Patient leakage current according IEC60601-1*		N.C. $I_{max} = 0,01$ mA S.F.C. $I_{max} = 0,05$ mA				
2.5 Equipment leakage current - alternative method according IEC62353:2007				$I_{max} = 1,00$ mA		
2.6 Applied part leakage current - alternative method according IEC62353:2007				$I_{max} = 0,05$ mA		
2.7 HF-leakage current bipolar according IEC60601-2-2:2009				$I_{max} = 55$ mA		
2.8 Insulation resistance power according IEC62353:2007 (Preparation according Service Manual)				test voltage = 500 V DC		
2.8.1 Mains against HF output				$R_{min} = 7$ M Ω		
2.8.2 Mains against housing				$R_{min} = 2$ M Ω		
2.8.3 Bipolar HF-output against housing				$R_{min} = 2$ M Ω		
2.9 HF-output power measurement with induction-free resistance						
HF-output		Operational mode	Dosage	R_{Load}	Rated value	
Bipolar		STANDART	60	75 Ω	60 Watt ± 20 %	
		FORCED	60	100 Ω	60 Watt ± 20 %	
2.10 DC resistance between the two HF output poles according IEC60601-2-2:2009 section 201.8.4.102				$R_{min} = 2$ M Ω		
2.11 Function test according to instructions for use					Performed:	
Place/Date		Inspector/Signatur		Responsible Organization		

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B BRAUN SHARING EXPERTISE		Inspection protocol – safety inspection				
TEST INTERVAL:		1 Year				
TYPE OF DEVICE:		Bipolar Electrosurgical Unit GN160				
MANUFACTURED BY:		Aesculap AG, Am Aesculap-Platz 78532 Tuttlingen/Germany				
SN:		INVENTORY NO.:		RESPONSIBLE ORGANIZATION:		
SCOPE OF INSPECTION					Finding	
1.) VISUAL INSPECTION					PASS	FAIL
1.1	Inspection of power cable for visually detectable damage					
1.2	Inspection of foot control for visually detectable damage					
1.3	Inspection of coagulation cable for visible signs of damage (record their number) and check the fitting of the plug connector					
1.4	Type plate still present and readable					
1.5	Check fuse links for prescribed values					
1.6	Overall condition of device: dirt, damage					
2.) ELEKTRICAL INSPECTION					Measured value	OK Yes / No
2.1	Protective earth resistance incl. mains cable according IEC62353:2007		$R_{max} = 0,3 \Omega$ (at $\geq 0,2 A$)			
2.2	Earth leakage current according IEC60601-1*		N.C.	$I_{max} = 0,50 mA$		
			S.F.C.	$I_{max} = 1,00 mA$		
2.3	Enclosure leakage current/touch current according IEC60601-1*		N.C.	$I_{max} = 0,10 mA$		
			S.F.C.	$I_{max} = 0,50 mA$		
2.4	Patient leakage current according IEC60601-1*		N.C.	$I_{max} = 0,01 mA$		
			S.F.C.	$I_{max} = 0,05 mA$		
2.5	Equipment leakage current - alternative method according IEC62353:2007		$I_{max} = 1,00 mA$			
2.6	Applied part leakage current - alternative method according IEC62353:2007		$I_{max} = 0,05 mA$			
2.7	HF-leakage current bipolar according IEC60601-2-2:2009		$I_{max} = 55 mA$			
2.8	Insulation resistance power according IEC62353:2007 (Preparation according Service Manual)		test voltage = 500 V DC			
	2.8.1	Mains against HF output	$R_{min} = 7 M\Omega$			
	2.8.2	Mains against housing	$R_{min} = 2 M\Omega$			
	2.8.3	Bipolar HF-output against housing	$R_{min} = 2 M\Omega$			
2.9	HF-output power measurement with induction-free resistance					
	HF-output	Operational mode	Dosage	R_{Load}	Rated value	
	Bipolar	STANDART	60	75 Ω	60 Watt $\pm 20 \%$	
		FORCED	60	100 Ω	60 Watt $\pm 20 \%$	
2.10	DC resistance between the two HF output poles according IEC60601-2-2:2009 section 201.8.4.102			$R_{min} = 2 M\Omega$		
2.11	Function test according to instructions for use				Performed:	
Place/Date		Inspector/Signatur		Responsible Organization		

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Inspection protocol – safety inspection

TEST INTERVAL: 1 Year
 TYPE OF DEVICE: Bipolar Electrosurgical Unit GN160
 MANUFACTURED BY: Aesculap AG, Am Aesculap-Platz 78532 Tuttlingen/Germany

RESponsible
 SN: INVENTORY NO.: ORGANIZATION:

SCOPE OF INSPECTION					Finding		
					PASS	FAIL	
1.) VISUAL INSPECTION							
1.1	Inspection of power cable for visually detectable damage						
1.2	Inspection of foot control for visually detectable damage						
1.3	Inspection of coagulation cable for visible signs of damage (record their number) and check the fitting of the plug connector						
1.4	Type plate still present and readable						
1.5	Check fuse links for prescribed values						
1.6	Overall condition of device: dirt, damage						
2.) ELEKTRICAL INSPECTION					Measured value	OK Yes / No	
2.1	Protective earth resistance incl. mains cable according IEC62353:2007		R _{max} = 0,3 Ω (at ≥0,2 A)				
2.2	Earth leakage current according IEC60601-1*		N.C. S.F.C.	I _{max} = 0,50 mA I _{max} = 1,00 mA			
2.3	Enclosure leakage current/touch current according IEC60601-1*		N.C. S.F.C.	I _{max} = 0,10 mA I _{max} = 0,50 mA			
2.4	Patient leakage current according IEC60601-1*		N.C. S.F.C.	I _{max} = 0,01 mA I _{max} = 0,05 mA			
2.5	Equipment leakage current - alternative method according IEC62353:2007		I _{max} = 1,00 mA				
2.6	Applied part leakage current - alternative method according IEC62353:2007		I _{max} = 0,05 mA				
2.7	HF-leakage current bipolar according IEC60601-2-2:2009		I _{max} = 55 mA				
2.8	Insulation resistance power according IEC62353:2007 (Preparation according Service Manual)		test voltage = 500 V DC				
	2.8.1	Mains against HF output		R _{min} = 7 MΩ			
	2.8.2	Mains against housing		R _{min} = 2 MΩ			
	2.8.3	Bipolar HF-output against housing		R _{min} = 2 MΩ			
2.9	HF-output power measurement with induction-free resistance						
	HF-output	Operational mode	Dosage	R _{load}	Rated value		
	Bipolar	STANDART	60	75 Ω	60 Watt ± 20 %		
		FORCED	60	100 Ω	60 Watt ± 20 %		
2.10	DC resistance between the two HF output poles according IEC60601-2-2:2009 section 201.8.4.102		R _{min} = 2 MΩ				
2.11	Function test according to instructions for use				Performed:		
Place/Date		Inspector/Signatur		Responsible Organization			

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Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany