DK50 D DK50 DM

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CE 0434

CONTENTS

1. GENERAL INFORMATION.....	2
PURPOSE.....	2
OPERATOR'S RESPONSIBILITY FOR PATIENT SAFETY	2
MARKINGS	2
WARNINGS.....	2
General safety warnings	3
Electrical system safety warnings	3
WARNING NOTICES AND SYMBOLS	4
USE.....	5
STORAGE AND TRANSPORT	5
2. EQUIPMENT DESCRIPTION.....	6
3. TECHNICAL DATA.....	8
4. OPERATION	9
INSTALLATION AND FIRST OPERATION	9
Preparing the compressor for use	9
Compressed air connection	10
Electrical connection	10
First operation	11
Accessories.....	11
OPERATION	12
Switching on the compressor	12
Running the compressor	13
Cleaning the suction filter.....	13
Cooling failure alarm	13
Cleaning the compressor	13
5. MAINTENANCE.....	13
REPAIRS AND SERVICE	13
MAINTENANCE SCHEDULE.....	15
NOTICE.....	15
Safety valve check	16
Replacing the suction filter	16
Replacing the blow down filter	16
Replacing the pressure regulator filter	17
Capacity check.....	17
Stabilizing the compressor before shipping.....	18
SHUTDOWN	18
EQUIPMENT DISPOSAL	18
6. TROUBLESHOOTING	19
7. SPARE PARTS	20
8. ELECTRIC AND PNEUMATIC DIAGRAMS	21
WIRING DIAGRAM	21
PNEUMATIC DIAGRAM	22
9. PARTS LIST	67

GENERAL INFORMATION

1. GENERAL INFORMATION

PURPOSE

The EKOM DK50 D, DK50 DM is a medical air compressor that supplies clean, oil-free compressed air for use with medical ventilators.

OPERATOR'S RESPONSIBILITY FOR PATIENT SAFETY

The Installation, Operation and Maintenance Manual is an integral part of the equipment and must be kept with the compressor. Careful review of this manual will provide information necessary for correct operation of the equipment.

Rx only

US Federal law restricts the sale of this device by or on the order of a physician.

MARKINGS

Products marked with the CE mark of compliance meet the safety requirements of the European Union (93/42/EEC).

WARNINGS

- The safety of operating personnel and trouble-free operation of the equipment are ensured only if original parts are used. Only accessories and spare parts mentioned in the technical documentation or expressly approved by the manufacturer can be used.
- If any other accessories or consumable materials are used, the manufacturer cannot be held responsible for the safe operation and functionality of the equipment.
- The warranty does not cover damages resulting from the use of accessories or consumable materials other than those recommended by the manufacturer.
- The manufacturer assumes responsibility for the safety, reliability and function of the equipment only if:
 - Installation, calibration, amendments, extensions and repairs are performed by the manufacturer, one of its representatives or a service provider authorized by the manufacturer
 - The equipment is used in accordance with the Installation, Operation and Maintenance Manual
- The Installation, Operation and Maintenance Manual accurately describes the design of the compressor and its compliance with safety and technical standards. The manufacturer reserves all rights to its wiring diagrams, procedures and names.

General safety warnings

The manufacturer developed and constructed the equipment so that damage would not occur when the equipment is used for its intended purpose. The manufacturer considers it its obligation to describe the following safety measures in order to avoid further damages.

- Equipment operation must comply with all local codes and regulations. In order to ensure safe work, the operator and user are responsible for following all applicable regulations.
- Original packaging should be kept for the possible return of the unit. Only original packaging ensures optimal protection of the equipment during transport. If it is necessary to return the equipment during the warranty period, the manufacturer is not liable for damages caused by incorrect packaging.
- The user must immediately notify the supplier if any problem occurs during the use of the equipment.
- This product is not intended for use in areas where there is a risk of an explosion. Do not operate the compressor in the presence of flammable anesthetics.
- Never feed oxygen or nitrous oxide into the compressor. Compressor components are not approved for oxygen or nitrous oxide use.
- This device can be used only with ventilators equipped with a low pressure alarm.

Electrical system safety warnings

- The equipment must only be connected to an appropriate power source that has correct grounding. In order to assure proper grounding, connect the compressor to a receptacle marked "hospital grade".
- Before the compressor is plugged in, make sure that the voltage and frequency of the mains specified on the equipment are the same as the power mains.
- Before operating, check for possible damage to the equipment and any connections. Damaged pneumatic and electrical lines must be replaced immediately.
- If a technical failure occurs, immediately disconnect the equipment from the mains (pull out the main power plug).
- During repairs and maintenance, ensure that:
 - The main power plug is removed from the power socket,
 - Compressed air lines are disconnected,
 - All pressure has been released from the air tank.
- Only a qualified technician can install this equipment.

GENERAL INFORMATION

WARNING NOTICES AND SYMBOLS

The following symbols are used for important information in the Installation, Operation and Maintenance Manual and on packaging and the product:

Attention, see instructions for use

Caution, risk of electric shock

Consult instructions for use

CE mark of compliance

Caution, hot surface

Handling mark on packaging – Fragile, handle with care

Handling mark on packaging – This way up (upright position)

Handling mark on packaging – Protect against moisture

Mark on packaging – Recyclable material

Handling mark on packaging – Temperature during storage and transport

Handling mark on packaging – Limited stacking


Alternating current

Ground connection

Equipotentiality

Fuse

USE

- The equipment can only be installed and operated in a dry, ventilated and dust-free area. Climatic conditions for operation - see Technical data.
- The compressor must stand on a flat and stable base.
- The compressor must not be exposed to rain. The equipment must not be used in humid or wet environments. Never use the compressor in the presence of flammable liquids or gases.
- Before connecting the compressor to respiration equipment, make sure that it meets the requirements of the respiration equipment. Refer to the Technical data for this purpose.
- Any use other than the compressor's intended use is not considered to be safe. The manufacturer is not responsible for any damages that result if the compressor is used for any other purpose. Risk is exclusively assumed by the operator/user.

STORAGE AND TRANSPORT

The compressor is shipped from the factory in transport packaging with the pump stabilized, protecting it from damage during transport.

For transport, always use the original packaging and secure the compressor in the upright position.

Protect the compressor from humidity, contamination and extreme temperatures during transport and storage. A compressor in its original packaging should be stored in a warm, dry and dust-free area.

Keep the packaging material, if possible. If not, dispose of the packaging material in an environmentally-friendly way. Cardboard can be recycled.

Before moving or transporting the compressor, the pressure in the air tank and hoses must be released and any condensed water must be drained. Secure the motor to prevent movement before shipping.

Prior to transport it is necessary to secure the motor inside the compressor (Chapter 5)


2. EQUIPMENT DESCRIPTION

1. PRESSURE - output air pressure gauge (accuracy $\pm 5\%$)
2. HUMIDITY - humidity indicator (performance)
3. Cooling failure indicator
4. Power indicator
5. Hour meter
6. Main input, power switch, primary fuses
7. Equipotential (ground) pin
8. Socket for the electrical cord
9. WALL - inlet air coming from an outside compressed air source
10. OUT - output air from compressor
11. Air tank
12. Suction filter
13. Pressure regulator with filter
14. Filter and water trap
15. Cooler
16. Control electronics
17. Safety valve
18. Internal back cover
19. Tube for supporting equipment

The equipment consists of an oil-free piston compressor driven by a single-phase electric motor, filtering and drying equipment, an air tank and a pressure regulator. The compressor is placed in a box lined with noise absorbing insulation. The equipment produces filtered compressed air without any traces of oil with reduced water content (dry air).

The air compressor has an input air connector labeled WALL (9) that allows the connection of an outside air source. When an outside compressed air source is connected to the WALL inlet fitting supplies air above the adjusted output pressure, the compressor will not be activated even though the power switch is in the ON position. Instead, the air from the outside source will be directed to the compressor's air output. When the outside air source pressure is less than the adjusted output pressure (or is not connected) and the power switch is in the ON position, the compressor is activated and supplies compressed air through the air output.

The compressor pulls in atmospheric air through the suction filter (12). Compressed air is cooled in the cooler (15), continues to the filter (14) and condensed liquid from the water trap is automatically released to the evaporation tray. The air continues through the check valve to the air tank (11) and through the pressure regulator with the output filter (13) to the OUT connector (10). Fans operate during compressor operation and when not operating when the temperature in the cabinet is over 40°C (104°F). If the temperature in the cabinet is over 80°C (176°F), a cooling failure is signaled by the red indicator (3) and an acoustic alarm. Control electronics (16) control compressor operation based on the pressure in the air tank (11). A safety valve (17) ensure the maximum allowed tank pressure is not exceeded.

TECHNICAL DATA

3. TECHNICAL DATA

TYPE	DK50 D					DK50 DM		
Output flow at pressure 3.5 bar (51 psig) Liters/min	40					60		
Output flow max. at pressure 3.5 bar (51 psig)	200 L/min for 2 s *							
			USA Model	100 /50-60			USA Model	
Rated voltage / Frequency V/Hz	230/50	230/60	120/60		230/50	230/60	120/60	
Nominal current A	2.8	2.8	5.6	5,6	3.4	4.3	8.6	
Filtration of air to µm	5							
Dew point depression at 40 L/min, 3 bar (43.5 psig), 20°C (68°F)	5°C (9°F) under the ambient temperature							
Outlet connection	RECTUS 21KS KS08 MPC Or DISS 1160-A (3/4"-16 UNF)							
Sound level dB(A)	≤ 49	≤ 51	≤ 51	≤ 51	≤ 50	≤ 52	≤ 52	
Mode of operation	Continuous - S1							
Release of condensed water	Automatic							
Indication of cooling failure	Acoustic and visual at 80°C (176°F)							
Output pressure	3.0 bar (43 psi) Adjustable to max. 3.5 bar (51 psi)							
Automatic start-up pressure (backup)	3.0 bar (43 psig) Adjustable to max. 3.5 bar (51 psi)							
Air tank capacity L	5							
Pressure range	5.0 - 7.0 bar (72.5 - 101.5 psi)							
Safety valve setting	8.0 bar (116 psi)							
Adjustment of the output air pressure	Pressure regulator							
Dimensions of compressor w x d x h	400x430x780 mm (16x17x31 in)				400x500x780 mm (16x20x31 in)			
Dimensions of packaging w x d x h	480x530x835 mm (19x20x33 in)				480x580x835 mm (19x23x33 in)			
Net weight	46 kg (101 lbs)				59 kg (130 lbs)			
Gross weight	60 kg (132 lbs)				73 kg (161 lbs)			
Implementation pursuant to EN 60601-1, EN 12021	Type B class I.							
Classification acc. to MDD 93/42 EEC, 2007/47 EC	II b							

* With 0.6 bar pressure drop (8.7 psi)

Climatic conditions for storage and transport

Temperature -25°C to +55°C (-13°F to +131°F), 24 hrs +70°C (+158°F)

Relative air humidity 10% to 90% (no condensation)

Relative air humidity for seaworthy packing 10% to 100% (with condensation)

Climatic conditions for operation

Temperature +5°C to +40°C (+41°F to +104°F)

Relative air humidity up to +95%

IPX0 Rating

4. OPERATION

INSTALLATION AND FIRST OPERATION

Do not use the compressor immediately after unpacking as it will not adjust to the ambient temperature.

Only qualified personnel can install the compressor and put it into operation for the first time. The installer shall train the operating personnel in the use and routine maintenance of the equipment. Installation and training of personnel should be acknowledged by the installer's signature on the installation certificate.

Prior to start-up, the four screws for transport stabilization must be removed. If the compressor is switched on without removing the stabilizers, the compressor could be permanently damaged!

The compressor does not contain a backup power supply.

This compressor can only be used with ventilators equipped with low pressure alarms.

Never block the suction filter and the vent outlets on the lower part of equipment. Never lift the unit by the handle! The handle is to be used to move a compressor on the floor with the brakes released.

If the compressor is equipped with a main source of air, the standby air source must be available.

ANY MODIFICATION OF THIS EQUIPMENT IS FORBIDDEN!


If this equipment is used nearby other instruments, the equipment must be observed in order to verify normal operations in the configuration it will be used.

Instruments may be affected electro-magnetically!

Preparing the compressor for use

After unpacking, place the compressor upright with the castor wheels firmly on the ground.

OPERATION

- Remove the front and back cover unscrewing the six screws on the cabinet side.
- Unscrew the four screws and remove the internal back cover (18).
- Remove all packing elements above (21) and under the compressor (20).
- Screw the internal back cover back on, connect the ground wire, screw the front and back cover and insert the suction filter.

Please retain the packing elements for future compressor transportation.

Compressed air connection

The medical compressor is equipped with quick couplings **WALL** (9) and **OUT** (10) to the rear of the cabinet.

Connect the pressure hose from the given equipment/respiration equipment to the quick coupling **OUT** (10) -outlet compressed air.

Connect the compressed air line from the central distribution line to the quick coupling **WALL** (9) -inlet compressed air. Air from the distribution system is automatically connected through the compressor to the OUT outlet port. In this configuration, the compressor serves as a backup source of compressed air. If the air pressure from the central distribution system is reduced, the compressor automatically switches on and there is no interruption in the supply of continuously pressurized air at the outlet of the compressor.

Please note that the air supply from the central distribution line entering the compressor must be medical grade air (particulate size, humidity.) The compressor does not modify the air from the central distribution line.

The hose connecting the compressor to the respiration equipment must not pass through a cold environment i.e. placed on the ground. It should be as short as possible with no kinks (this may cause water to condense inside the hose).

Electrical connection

The compressor comes with a plug containing an appropriate protective contact (ground.) Adhere to local electrical regulations. The voltage and frequency of the mains must comply with the specifications on the data label.

The electrical cord must not be stressed or have any tension exerted upon it, and must always be free.

- The socket must be accessible for safety reasons so that the equipment can be safely disconnected in case of an emergency.
- The relevant current circuit must be protected.

- Connection of the ground connection (7) to other equipment must adhere to local electrical regulations.
- Fasten the electrical cord through the holder (8).

First operation

- Make sure the stabilizing screws used during transport were removed.
- Check that the connection to the compressed air supply is correct.
- Check for proper connection to the main power supply.
- Switch on the pressure switch (6) to position "I".

After the compressor is first put into operation, the air tank of the compressor is pressurized to a 'switching-off' pressure and the compressor will then be automatically switched off. During the following time period the compressor works in an automatic mode, i.e. according to the consumption of compressed air; the pressure switch controls compressor operation. During operation, the equipment releases condensed water from the condenser (filter) automatically to the evaporation tray.

Accessories

Holder for supporting equipment - Item 604011346

The holder (19) is suitable for supporting a respiratory apparatus with accessories. Please remove the plug in the middle of the upper side before installation. Screw the holder into the socket and tighten it fast.

Technical data:

Diameter 25 mm
Height 1100 mm above the floor

The maximum holder load is 28 kg. The centre of gravity of supporting equipment must be in the axis of the holder, max. 1100 mm above the floor. Consult the manufacturer regarding any other connections.

The manufacturer is not liable for any damage resulting from excessive loads on the equipment!

The supplier is obliged to ensure the acceptability of a load delivered by a supplied accessory.

It is forbidden to lean or press against an installed accessory for support!


The equipment must be lifted when travelling over an obstacle!
Supporting equipment must be disassemble before transport!

OPERATION

Kit for higher relative air humidity - Item 604011497

1. Hose	PUR 6/4-130	062000337	1
2. Hose	PUR 6/4-700	062000357	1
3. Tank	LUH-1	062000333	1
4. Holder	3KB-972	023000564	1
5. Washer	D5 043000002, 043000005		4
6. Screw	M5 041000131		4
7. Reducer	10/6 062000339		1
8. Connector	L6/6 025400153, 025400154		1
9. Clamp	self-adhesive 033400011		2

Dismount the hose (A) and plug (B). Install the reducer (7), connector (8) and plug-in hose (1). Screw holder (4), plug-in hose (2) and install tank (3). Fix hose with two clamps (9).

Empty tank regularly.

OPERATION

In case of emergency, switch the equipment off at the switch and pull out the main power plug.

Automatic operation. When the pressure in air tank falls to the 'switch on' pressure, the compressor is automatically switched on. The compressor is automatically switched off when the pressure increases to the 'switch off' pressure in the air tank.

Switching on the compressor

The compressor is switched on at the main power switch (6) by putting it in position "I". Switching on is signaled by the green indicator (4).

Running the compressor

When the pressure of the air at the 'Inlet' port is reduced or if there is a reduction of pressure in the air tank, the compressor is switched on and the air tank is filled with compressed air. After the 'switch-off' pressure is reached, the compressor is automatically switched off. Pressure of the output air from compressor is indicated on the pressure gauge (1).

During compressor operation, the pointer of the humidity indicator (2) should be in the green zone. If equipment air consumption is higher than maximum output, the pointer of indicator will move from the green zone and the air from compressor will not reach the declared parameters (see Technical data).

The hour meter (5) indicates total running time of the compressor.

Cleaning the suction filter

At least once a week take out and clean the suction filter (12) located on the front side. Wash the filter in warm soapy water, rinse thoroughly and allow it to dry. Insert the clean filters so that the suction openings are completely covered by the filters.

Cooling failure alarm

Problems with compressor cooling are signaled by the red indicator (3) and the acoustic alarm.

If a cooling failure alarm is signaled, secure an alternate air pressure source and disconnect the compressor.

Cleaning the compressor

The product can be cleaned using normal detergents. To clean the compressor, use a detergent that contains no abrasives, chemical solvents or other corrosive agents.

5. MAINTENANCE

REPAIRS AND SERVICE

Warranty and extended warranty repairs are to be completed by the manufacturer or a service provider authorized by the manufacturer.

The manufacturer reserves the right to modify the equipment in any way that will not alter the function or the operation of the equipment.

OPERATION

Only a qualified technician or the Customer Service Department of the manufacturer may perform repairs that go beyond routine maintenance. Use only spare parts and accessories approved by the manufacturer.

Prior to any maintenance or repairs, switch off the compressor and disconnect it from the mains (pull out the main power plug).

MAINTENANCE SCHEDULE**NOTICE**

The operating entity is obliged to ensure that all tests of the equipment are carried out repeatedly at least once within every 24 months (EN 62353) or in intervals as specified by the applicable national legal regulations. A report must be prepared on the results of the tests (e.g.: according to EN 62353, Annex G), including the measurement methods used.

Maintenance range	Page	Time interval	To be performed by
Clean suction filter	12	At least once a week	Staff
Safety valve check	14	Once per year	Qualified expert
Replace filter (B)	14	Every 8000 hours	Qualified expert
Replace filter in blow down filter (C)	14	Every 2000 hours	Qualified expert
Replace filter in the pressure regulator (D)	15	Every 2000 hours	Qualified expert
Test the tightness of joints and inspect the equipment	Service documentation	After two years	Qualified expert
Replace piston complete with piston rings, O-rings and bearing	Service documentation	Every 4000 hours	Qualified expert
Bearing replacement	Service documentation	Every 4000 hours	Qualified expert
Perform "Repeated Test" according to EN 62353	13	1 x 2 years	Qualified expert


MAINTENANCE

Safety valve check

The safety valve is adjusted to 8 bar (116 psig) by the manufacturer, then tested and marked. It must not be readjusted.

- Turn the screw of the safety valve (17) several rotations to left until the safety valve puffs.
- Only allow the safety valve freely puff **for a short period of time**.
- Turn the screw to the right until it is tight, closing the valve.


Never use the safety valve to release pressure in the air tank. This could damage the safety valve.

Replacing the suction filter

- Unlock the cover and pull it out.
- Replace the filter
- Lock the cover

Spare part:

Suction filter 05W POLYESTER, Item 025200194

Replacing the blow down filter


- Loosen lock on filter (14) by pulling it down.
- Slightly turn the cover of the filter and pull it out.
- Unscrew the filter holder.
- Replace filter (C), screw in the filter holder.
- Place the cover of filter and secure it by turning to right till it locks.

Replacing the pressure regulator filter

Before working on this equipment all air pressure must be released from the tank so that pressure is zero and the equipment must be disconnected from the electrical power mains.

- Unscrew the filter regulator cover (13) and pull it out.
- Unscrew the filter holder.
- Replace filter (D), screw in the filter holder.
- Screw in the regulator cover.

Setting output air pressure

Switch on the compressor and set the output pressure air flow until the humidity indicator (2) remains in the green area. Unlock the control button of the filter regulator (13) and rotate to set the demand output pressure. After setting the pressure, lock the control button of the regulator by pushing down.

Capacity check

Before checking capacity it is necessary to switch off the compressor using the power switch, vent air pressure out of the air tank and disconnect the output hose from the OUT (10) fitting. Next switch the compressor on using the power switch and measure (record) the time between start and stop (2). Compressor capacity is met when the recorded time is not more than 45 seconds for DK50 D or 30 seconds for DK50 DM.

MAINTENANCE

Stabilizing the compressor before shipping

Before shipping, the compressor must be stabilized to prevent movement. The procedure is illustrated on page 9.

- Remove the front and back covers unscrewing screws on the cabinet side, disconnect the ground connector on the covers.
- Disassemble the internal back cover (18) unscrewing four screws.
- Insert the shipping elements (21) under and above the compressor.
- Assemble the internal back cover (18), connect the ground connector on the covers, and insert the covers and the suction filter.

SHUTDOWN

If the compressor is not going to be used for a long period of time, disconnected it from the main power supply and release the pressure in the air tank.

EQUIPMENT DISPOSAL

- Disconnect the equipment from the main power supply.
- Release the pressure in the air tank.
- Dispose of the equipment according to local regulations.
- Parts used in this product have no negative impact on the environment when disposed of properly.

6. TROUBLESHOOTING

Prior to servicing the equipment, reduce the pressure in the air tank to zero and disconnect the equipment from the main power supply.

Only trained service personnel can perform the activities listed in the troubleshooting guide.

PROBLEM	POTENTIAL CAUSE	TROUBLESHOOTING
Compressor does not start	No main power voltage	Switched off main breaker in distribution system
	No current from the pressure switch	Check supply
		Replace a defective fuse
		Loose wire terminal – tighten
		Power cord defective – replace
	Interrupted winding of motor, damaged thermal protection	Replace motor
	Defective capacitor	Replace capacitor
The humidity indicator is out of the green zone	Jammed piston or another parts	Replace damaged parts
	Pressure switch does not switch on	Change pressure switch
Compressor often switches on	Leakage of air from pneumatic distribution system	Inspect pneumatic distribution system – seal loose connection
	Large consumption of compressed air	Do not exceed max. flow (see Technical data)
The operation of compressor is extended	Leakage of air from pneumatic distribution system	inspect pneumatic distribution system – seal loose connection
	Leakage of check valve	Clean check valve, replace seal, replace check valve
Red indication lights – state of overheating	Leakage of air from pneumatic distribution system	Inspect pneumatic distribution system – seal loose connection
	Worn out piston rings	Replace worn out piston rings
	Dirty filter inserts	Replace the dirty inserts
Compressor is noisy (knocking, metal noises)	Dirty suction filter	Clean or replace the dirty suction filter
	Dirty filter inserts	Replace the dirty inserts with new ones
	Cooling fans don't work	Replace fans or faulty thermal switch
	Unit is hot, unventilated area	Relocate unit
Water coming out of outlets	Damaged piston bearing or bearing on motor	Replace damaged bearing
	Loose (cracked) belt on the compressor hanger	Replace damaged hanger
Water coming out of outlets	Water trap (13) or (14) not working	Clean or replace the water trap

7. SPARE PARTS


• Suction filter (12)	4KB-552	025000011
• Suction filter (B)	05W POLYESTER	025200194
• Filter (C)	111511-40B	025200070
• Filter (D)	AF20P-060S	025200113
• Fuse		
• DK50 D 230V	T6.3A	038100004
• DK50 D 120V, DK50 DM 230V	T10A	
	038100005	
• DK50 DM 120V	T16A	038100007
• Nipple		
• RECTUS	21SSTF08MPC	025500131
• DISS	1160-A	024000261

8. ELECTRIC AND PNEUMATIC DIAGRAMS

WIRING DIAGRAM


1/N/PE ~ 115/230 V 50/60 Hz
ELECTRIC OBJECT OF 1st CAT.
B TYPE

DK50 D

1/N/PE ~ 115/230 V 50/60 Hz
ELECTRIC OBJECT OF 1st CAT.
B TYPE

DK50 DM


Note

Wire entering the terminal should be marked with marking element with corresponding number.

Des.	Name	Des.	Name
HL1,HL2	Glow lamp	M	Electric motor
EV1,EV2	Ventilator	ST1,ST2	Thermal switch
PH	Hour meter	XC	Socket with fuses
YV	Solenoid valve	Cb1	Capacitor

ELECTRIC AND PNEUMATIC DIAGRAMS

PNEUMATIC DIAGRAM

1. Suction filter
2. Suction filter
3. Compressor
4. Solenoid valve
5. Sound absorber
6. Cooler
7. Micro filter with water trap
8. Check valve
9. Air tank
10. Pressure regulator with filter
11. Air Outlet
12. Air Inlet
13. Main input, power switch, primary fuses
14. Pressure sensor
15. Humidity indicator
16. Pressure gauge
17. Pressure relief valve
18. Fan
19. Water tray
20. Shuttle valve

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