BSC Detachable Coil Power Supply

| | |
|------------------------------------|----|
| Directions for Use | 3 |
| Instrucciones de uso | 11 |
| Mode d'emploi | 19 |
| Gebrauchsanweisung | 27 |
| Istruzioni per l'uso | 35 |
| Gebruiksaanwijzing | 43 |
| Instruções de Utilização | 51 |

BSC Detachable Coil Power Supply

Directions for Use

CAUTION
Federal Law (USA) restricts this device to sale by
or on the order of a physician.

I. INTRODUCTION / DEVICE DESCRIPTION

The BSC Detachable Coil Power Supply is a battery-operated instrument designed to detach Boston Scientific's GDC™ and Matrix™ Detachable Coils. GDC and Matrix™ Detachable Coils are connected to the distal end of a stainless steel delivery wire.

Note: The BSC Detachable Coil Power Supply is not intended for use with Target™ Detachable Coils.

Using fluoroscopic guidance, the physician inserts a Boston Scientific 2-tip marker Infusion Catheter into the femoral artery of the patient and guides the distal tip to the aneurysm. A GDC or Matrix™ Detachable Coil is then inserted into the catheter, and the physician guides the coil to the desired position inside the aneurysm before detachment is initiated.

The black (negative) terminal of the BSC Detachable Coil Power Supply is connected to a return electrode on the patient; a sterile 20 or 22 gauge stainless-steel hypodermic uncoated needle inserted at the groin. The red (positive) terminal of the BSC Detachable Coil Power Supply is connected to the proximal end of the delivery wire.

The BSC Detachable Coil Power Supply is then turned on. After passing self-diagnostic tests, a DC current of 1 mA is sent through the delivery wire and returns through the patient. The current may be adjusted to a different setting by the physician at any time by pressing the current select switch  on the front panel.

The current flowing through the delivery wire will begin to dissolve, by electrolysis, the exposed stainless steel at the junction with the platinum coil. Over a period of several minutes, the exposed steel will completely dissolve and the platinum coil detaches. When the coil detachment is detected, the BSC Detachable Coil Power Supply immediately halts current flow, freezes the displays and emits five beeps to indicate to the physician that coil detachment has occurred.

Using fluoroscopy, the physician can verify that the delivery wire may be withdrawn without disturbing the newly placed coil. If further electrolysis is necessary, the procedure can be resumed by pressing the current select switch .

II. INTENDED USE / INDICATIONS FOR USE

The Boston Scientific Corporation (BSC) Detachable Coil Power Supply is intended for use with Boston Scientific's GDC Detachable Coils and Matrix™ Detachable Coils in the embolization of intracranial aneurysms and other vascular malformations of the neuro and peripheral vasculature.

Refer to GDC or Matrix™ Detachable Coil
Directions For Use for use in the clinical setting.

III. CONTRAINDICATIONS

Refer to GDC or Matrix™ Detachable Coil Directions For Use.

IV. WARNINGS

Refer to GDC or Matrix™ Detachable Coil Directions For Use.

This device should only be used by physicians who have received appropriate training in interventional neuroradiology or interventional radiology and preclinical training on the use of this device as established by Boston Scientific.

The use of cables and/or other medical devices other than the BSC Detachable Coil Connecting Cables specified may result in increased emissions or decreased immunity of the BSC Detachable Coil Power Supply.

The BSC Detachable Coil Power Supply should not be used adjacent to or stacked with other equipment; if adjacent or stacked use is necessary, the BSC Detachable Coil Power Supply should be observed to verify normal operation in the configuration in which it will be used.

V. CAUTIONS

A second, back-up BSC Detachable Coil Power Supply is required for all procedures.

The BSC Detachable Coil Power Supply needs special precautions regarding EMC. Use of the BSC Detachable Coil Power Supply should be according to the EMC information provided within these Directions for Use.

Portable and mobile RF communications equipment may affect the performance of the BSC Detachable Coil Power Supply.

VI. ADVERSE EVENTS

Please notify your Boston Scientific representative immediately by telephone or FAX if a device malfunctions or patient complication or injury is experienced or suspected associated with the use of this device. Please make every attempt to retain any suspect device, its associated components and their packaging for return to Boston Scientific.

No adverse events are associated with the use of the BSC Detachable Coil Power Supply as a stand alone device. Refer to GDC™ or Matrix™ Detachable Coil Directions for Use for adverse events associated with the use of the BSC Detachable Coil Power Supply with GDC or Matrix Detachable Coils.

VII. CONFORMANCE TO STANDARDS

EMF output conforms with EN 60601-1-2:2007, Medical electrical equipment - Part 1-2: "General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests."

VIII. TECHNICAL DESCRIPTION

BSC DETACHABLE COIL POWER SUPPLY SPECIFICATIONS

The BSC Detachable Coil Power Supply operates with the following:

| | |
|--|---|
| Power | 9 VDC from two (2) parallel 9V batteries |
| Batteries | 9V Alkaline Batteries: Two (2) Eveready, Part No. 522 or Two (2) Duracell, Part No. 1604B2 |
| Expected Battery Life | 5 hours |
| Red Cable | Positive (anode) |
| Black Cable | Negative (cathode) |
| Current | 1 mA Max. |
| Voltage | 12 VDC Max. |
| Operating Temp. | 5° to 45° C (41° to 113° F) |
| Storage Temp. | -20° to 60° C (-4° to 140° F) |
| Relative Humidity | 10% to 90% (non-condensing) |
| Unit Size | 23.50 cm x 12.07 cm x 6.35 cm (9.25" x 4.75" x 2.5") |
| Unit Weight | 0.91 kg (2 lbs) |
| Power Source Location (EN 60601-1) | Internally Powered Equipment |
| Device Classification (EN 60601-1) | Type BF |
| Compatibility with Flammable Anaesthetic Mixtures (EN 60601-1) | Equipment not suitable for use in the presence of a Flammable Anaesthetic Mixture with Air or with Oxygen or Nitrous Oxide |
| Protection against Ingress of Water (EN 60601-1) | IPX0 |
| Mode of Operation (EN 60601-1) | Continuous Operation |

Guidance and manufacturer's declaration - electromagnetic emissions

The BSC Detachable Coil Power Supply is intended for use in the electromagnetic environment specified below. The customer or the user of the BSC Detachable Coil Power Supply should assure that it is used in such an environment.

| Emissions test | Compliance | Electromagnetic environment - guidance |
|--|----------------|--|
| RF emissions CISPR 11 | Group 1 | The BSC Detachable Coil Power Supply uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF emissions CISPR 11 | Class B | The BSC Detachable Coil Power Supply is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2 | Not applicable | |
| Voltage fluctuations/flicker emissions IEC 61000-3-3 | Not applicable | |

| Guidance and manufacturer's declaration - electromagnetic immunity | | | |
|---|---|--|---|
| The BSC Detachable Coil Power Supply is intended for use in the electromagnetic environment specified below. The customer or the user of the BSC Detachable Coil Power Supply should assure that it is used in such an environment. | | | |
| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment - guidance |
| Electrostatic discharge (ESD) IEC 61000-4-2 | ± 6 kV contact ± 8 kV air | ± 2, 4, 6 kV contact ± 2, 4, 8 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. |
| Electrical fast transient/burst IEC 61000-4-4 | ± 2 kV for power supply lines ± 1 kV for input/output lines | Not applicable | Not applicable |
| Surge IEC 61000-4-5 | ± 1 kV line(s) to line(s) ± 2 kV line(s) to earth | Not applicable | |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | <5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 s | Not applicable | |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 3 A/m | Not applicable | |
| NOTE: U_T is the a.c. mains voltage prior to application of the test level. | | | |

| Guidance and manufacturer's declaration - electromagnetic immunity | | | |
|---|---|-----------------------------|---|
| The BSC Detachable Coil Power Supply is intended for use in the electromagnetic environment specified below. The customer or the user of the BSC Detachable Coil Power Supply should assure that it is used in such an environment. | | | |
| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment - guidance |
| Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3 | 3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz | Not applicable 3 V/m | Portable and mobile RF communications equipment should be used no closer to any part of the BSC Detachable Coil Power Supply, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance Not applicable $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:  |
| NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. | | | |
| NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. | | | |
| ^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the BSC Detachable Coil Power Supply is used exceeds the applicable RF compliance level above, the BSC Detachable Coil Power Supply should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the BSC Detachable Coil Power Supply. | | | |

| Recommended separation distances between portable and mobile RF communications equipment and the BSC Detachable Coil Power Supply | | | |
|--|--|--|---|
| The BSC Detachable Coil Power Supply is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the BSC Detachable Coil Power Supply can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BSC Detachable Coil Power Supply as recommended below, according to the maximum output power of the communications equipment. | | | |
| Rated maximum output power of transmitter W | Separation distance according to frequency of transmitter m | | |
| | 150 kHz to 80 MHz | 80 MHz to 800 MHz $d = 1.2\sqrt{P}$ | 800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$ |
| 0.01 | Not Applicable | 0.12 | 0.23 |
| 0.1 | Not Applicable | 0.38 | 0.73 |
| 1 | Not Applicable | 1.2 | 2.3 |
| 10 | Not Applicable | 3.8 | 7.3 |
| 100 | Not Applicable | 12 | 23 |
| For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. | | | |
| NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. | | | |
| NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. | | | |

IX. TECHNICAL DRAWING


Figure 1
BSC Detachable Coil Power Supply Controls & Displays

X. CONTROLS & DISPLAYS

POWER ON/OFF BUTTON

This momentary-action button is used to turn the BSC Detachable Coil Power Supply on and off. Pressing the button one time turns the unit on, and the unit runs a self-diagnostic check for three seconds. Once completed, the BSC Detachable Coil Power Supply takes 10-15 seconds to reach the selected output current setting.


Pressing the button again will turn the unit off instantly. **Note:** In Pause mode, the unit will automatically turn off after 15 minutes have elapsed.

CURRENT SETTING/RESUME CONTROL

This momentary-action, touch-sensitive membrane switch is used to change the patient current setting. Each time the BSC Detachable Coil Power Supply is turned on, the current is set to 1.00 mA. Pressing the membrane one time changes the setting to 0.5 mA; pressing a second time changes the setting to 0.75 mA; pressing a third time returns the setting to 1.00 mA. Each time the switch is pressed, the current display briefly flashes the new current setting.

The *resume control* is designed to restart the flow of current through the system without resetting the elapsed detachment time. For instance, if the BSC Detachable Coil Power Supply incorrectly signals coil detachment, current flow can be reestablished by simply pressing the current select switch.


LED DISPLAYS & INDICATORS

- **Voltage LED Display & Indicator** 

The three-digit, red LED voltage display is located on the top of the front panel and displays the output voltage with an accuracy of 0.01 volts DC.
- **Current LED Display** 

The three-digit, red LED current display is located in the center of the front panel and displays the current flowing to the GDC™ or Matrix™ Detachable Coil with an accuracy of 0.01 mA. In addition, the display briefly flashes the new current setting (such as "1.00" mA) when the front panel current select switch is pressed. The readout then returns to the continuous display of actual current.
- **Time LED Display** 

The four-digit, red LED time display is located at the bottom of the front panel and displays in minutes and seconds. The time readout is the elapsed time that current has been flowing through the BSC system.
- **Check Indicator** 

The CHECK indicator light will flash yellow ON and OFF at the beginning of a detachment procedure if there is a poor connection to the patient. It is located in the lower left corner of the front panel.
- **Detach Indicator** 

The DETACH indicator light flashes ON and OFF when the unit has detected detachment of the coil from the delivery wire. This indication is accompanied by five beeps. It is located on the bottom center of the front panel.
- **Battery Indicator** 

The BATTERY indicator light illuminates red when the unit is turned on and when the voltage of the internal 9-volt batteries drops below approximately 7.5V. It is located on the lower right corner of the front panel. New batteries should be installed prior to each treatment session.

PAUSE MODE

The Pause mode is designed to allow the clinician, utilizing fluoroscopy, to confirm coil detachment prior to turning off the BSC Detachable Coil Power Supply. Once the BSC Detachable Coil Power Supply circuitry detects coil detachment, the unit will emit five (5) audible beeps signaling detachment. In addition, current flow stops and all displays freeze, showing indicated levels at the time of detachment. At this point, the unit enters the Pause mode.

XI. DIRECTIONS FOR USE / DETACHMENT PROCEDURE

Refer to GDC or Matrix[®] Detachable Coil
Directions For Use for use in the clinical setting.

XII. TROUBLESHOOTING HINTS

- **CHECK**  **indicator flashes when unit is first turned on.**

If the CHECK indicator flashes ON and OFF when the unit is first turned on, a poor connection most likely exists at the site of the return electrode (a sterile 20 or 22 gauge stainless-steel hypodermic uncoated needle inserted at the groin) or at the positive (red) cable connection to the GDC or Matrix[®] Detachable Coil delivery wire. Stop the procedure, depress the On/Off button, and reconnect the connecting cables. If the problem persists, replace the connecting cables. If the problem continues to persist, use the back-up BSC Detachable Coil Power Supply and begin the procedure again.
- **Current**  **delivery ceases for 2 seconds and resumes at selected current level.**

Occasionally, the voltage required to maintain a constant, controlled current increases significantly due to changes in system impedance without an actual coil detachment. When this occurs, the system will initiate a "rest cycle" in an attempt to reduce the system impedance. This "rest cycle" will only be initiated once during a detachment procedure. If a second, similar increase in voltage is observed, the system will signal a detachment, which must be verified under fluoroscopy.
- **LED readouts partially visible.**

If any of the displayed numbers are misformed, stop the procedure, and depress the On/Off button. During initial power-on, the displays should first appear blank and subsequently display a bright red "8." If the problem persists, restart the procedure by using the back-up BSC Detachable Coil Power Supply and notify Boston Scientific immediately.
- **ERROR MESSAGE: "CPU. FAIL"**

This message is displayed if one of the following tests fails during the power-on diagnostics:

 - a. the Program Integrity (EEPROM checksum) test, or
 - b. the internal RAM test.
- **ERROR MESSAGE: "OP. FAIL"**

This message is displayed if one of the following tests fails during the power-on diagnostics:

 - a. the auto-shutdown current-enable hardware test has failed, or
 - b. the output voltage is out of range.

This message is also displayed if the output voltage rises above 11.7 volts during a procedure.

- **ERROR MESSAGE: “Cur. FAIL”**

This message is displayed if the following test fails during the power-on diagnostics:

- a. the output current is out of range.

This message is also displayed if the output current is out of range during the procedure.

- **ERROR MESSAGE: “PIC. FAIL”**

This message is displayed if the following test fails during the power-on diagnostics:

- a. the patient isolation relay has failed.

CLEARING ERRORS

To clear an error, the BSC Detachable Coil Power Supply must be turned off and then on again. If the above error messages persist after restart, stop the procedure and disconnect the BSC Detachable Coil Power Supply. Continue the procedure by using the back-up BSC Detachable Coil Power Supply, and notify Boston Scientific immediately.

XIII. MISCELLANEOUS

REQUIRED ADDITIONAL ITEMS & ACCESSORIES

In addition to the GDC™ or Matrix™ Detachable Coils and two BSC Detachable Coil Power Supply systems (one for back-up), the following are required for each procedure:

- Two fresh batteries
- BSC Detachable Coil Connecting Cables [use only BSC Detachable Coil Connecting Cable sets which have a red cable with a shrouded-type connector pin (catalog number 451102-4)] – sold separately
- Return electrode (a sterile 20 or 22 gauge stainless-steel hypodermic uncoated needle inserted at the groin)
- Alcohol-dampened gauze to clean the proximal end of the GDC or Matrix® Detachable Coil delivery wire prior to attaching the cable

MOUNTING

Pole

The BSC Detachable Coil Power Supply can be mounted on a 1.9 - 3.8 cm (0.75" - 1.5") diameter pole by pulling back on the rear clip, aligning the unit's cylindrical cavity with the pole and gently releasing. **Note:** if IV pole is too narrow, an attachment adapter is packaged with the BSC Detachable Coil Power Supply.

Tabletop

The BSC Detachable Coil Power Supply unit can stand alone by pulling out the swing arm at the back of the unit.


Figure 2
Table & Pole
Mounting Positions

IV POLE CLAMP

Instructions for use:

1. Mount the two sections of the IV clamp (items 1 and 2 in Figure 3) around the IV pole.
2. Secure the sections together using the knurled screws (item 3).
3. Tighten screws securely.
4. Mount the BSC Detachable Coil Power Supply above the clamp and make sure the guiding shoulders of the two units interlock.


Figure 3
IV Pole Clamp

BATTERY CHANGE & DISPOSAL

Remove the battery access door. Replace used batteries with two new 9V alkaline batteries in the compartment as shown, noting **CORRECT POLARITY**. Snap battery door securely into place. Dispose of used batteries in accordance with hospital, administrative and/or local government policy. Remove batteries when unit is not in use.

Precaution: Incorrect placement of batteries will prevent unit from operating.

WHEN SYSTEM NOT IN USE

If either the BSC Detachable Coil Power Supply or the backup BSC Detachable Coil Power Supply is not used in any six-month period, perform the preventive maintenance procedure outlined below.


Figure 4
Battery Compartment for BSC Detachable
Coil Power Supply

Preventive Maintenance should be performed
twice per year if otherwise not in use.

Attach a 2.2k Ω , 5% tolerance resistor which has been placed in parallel with a 0.0022 μ F, 5% tolerance capacitor across the two jacks. Turn on unit at the default 1mA level and allow the self-check systems to run (approximately 10 seconds). Unit should read 2.2 V \pm 0.2 V. Leave system on for at least 10 seconds. To simulate coil detachment, break circuit by simultaneously disconnecting the resistor and the capacitor from the jacks. Unit should beep 5 times and the DETACH indicator should begin flashing ON and OFF. If unit does not function properly, notify Boston Scientific immediately.

HOW SUPPLIED

Do not use if package is opened or damaged

Do not use if labeling is incomplete or illegible.

Handling and Storage:

Store in a cool, dry, dark, place.

Each Boston Scientific carton contains one BSC Detachable Coil Power Supply, one IV stand clip-on adapter, and one Directions for Use.

After opening, inspect contents. Do not use the BSC Detachable Coil Power Supply if there are missing items or if there is evidence of damage. Report any problems immediately to Boston Scientific.

Do not destroy or discard original packing materials or container. Save for storage or re-shipping.

CLEANING

CAUTION

Do not sterilize BSC Detachable Coil Power Supply.

CAUTION

**Do not use water to clean BSC Detachable Coil Power Supply.
Do not submerge BSC Detachable Coil Power Supply.
Damage to the system could result.**

Use only a cloth dampened with alcohol to clean the plastic surface of the BSC Detachable Coil Power Supply. Avoid all types of detergents, cleaners, or solvents which might contain chemicals that could seriously damage the plastic or leave contaminating residues.

SERVICE

Service, if necessary, must be performed by qualified personnel at Boston Scientific. If you experience difficulty with the BSC Detachable Coil Power Supply, you may return it to Boston Scientific Customer Service Department for repair or replacement. Boston Scientific will, at its option, repair or replace the faulty unit. To return a BSC Detachable Coil Power Supply, obtain a Returned Goods Authorization number from Boston Scientific's Customer Service Department. Return the BSC Detachable Coil Power Supply with a written explanation of the reason for return. Remove batteries from the BSC Detachable Coil Power Supply, then pack it carefully in its original carton or any sturdy carton with enough packing material to prevent damage. Be sure to insure your package for protection.

CAUTION

**Do not remove cover (or back) from BSC Detachable Coil Power Supply.
Refer servicing to qualified Boston Scientific service personnel.**

If you have further questions, please contact Boston Scientific.

DISPOSAL

Dispose of cables and BSC Detachable Coil Power Supply in accordance with hospital, administrative and/or local government policy.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

PRECAUTION

The safety and performance characteristics of the BSC System (BSC Detachable Coils, BSC Detachable Coil Power Supply, delivery systems and accessories), when used with another manufacturer's devices (whether coils, coil delivery devices, catheters, guidewires, and/or other accessories) has NOT been established. Due to the potential incompatibility of non-Boston Scientific components with the BSC System, the use of another manufacturer's device(s) with the BSC System is not recommended.

Symbol Translation Key

| | |
|---|--|
|  | Dangerous Voltage |
|  | See Manual |
|  | Type BF Applied Part |
|  | cTUVus Mark indicates compliance to UL 60601-1 and CAN/CSA 22.2 601.1 M90 covering electrical safety requirements for the US and Canada. |
|  | Contents |
|  | Separate Collection |
| UPN | Product number |
|  | Date of Manufacture |
| SN | Serial Number |
| REF | Catalog Number |
|  | Consult instructions for use |
| AUS | Australian Sponsor Address |
|  | Legal Manufacturer |
| EC REP | EU Authorized Representative |
|  | Do Not use if package is damaged. |
|  | Recyclable Package |

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© 2011 Boston Scientific Corporation or its affiliates. All rights reserved.
Boston Scientific International S.A.
Boston Scientific (Australia) Pty Ltd