

More comfortable mammograms. More satisfied patients.

Hologic addresses the main reason too many women put off having a potentially life-saving mammogram: discomfort.¹ Our SmartCurve™ breast stabilization system delivers an exceptional patient experience without compromising image quality, workflow or dose.



The SmartCurve™ system is clinically proven to deliver a more comfortable mammogram without increasing dose.²



Curved compression surface mirrors the shape of a woman's breast.



95% of women surveyed would recommend facilities that use the SmartCurve™ system.²



The SmartCurve™ system will change the way women look at breast screening.

The fear of pain stops many women from making regular breast screening a healthcare priority. Not only does this put them at risk of a cancer being missed or diagnosed at a more advanced stage, it can also have a negative impact on patient volume¹ and satisfaction scores. The SmartCurve™ system was designed to address this challenge. It transforms the breast screening experience – for you and your patients.



SmartCurve™ compression surface conforms to the shape of the breast and is curved at the chest wall and sides to reduce pinching.

Product information

The SmartCurve™ system is intended for standard screening views; paddles will come in two sizes to accommodate a majority of women and breast sizes. It may not be appropriate for all procedures. Refer to the Dimensions Platform datasheet for guidelines on usage.

Ordering details

| Part Number | Description |
|----------------|--|
| DIM-ASY-SCS-V2 | SmartCurve™ Breast Stabilization system for Selenia® Dimensions® systems |

References

1. Padoan M, Ferrante D, Pretti G, Magnani C. Study of socio-economic characteristics, diagnosis and outcome of women participating or not participating in mammogram screening. *Ann Ig.* 2014;26(6):518–526. 2. Data On File: DHM-06039 Rev 002 Jennifer Bartoshevich, Teri Orefice, Sajjad Mansoor, et.al Internal Study comparing Hologic's flat paddle to the SmartCurve paddle (18x24cm) (2017).

SS-00523-EUR-EN Rev.003 (02/20) Hologic Inc. ©2020 All rights reserved. Hologic, Dimensions, MammoPad, Selenia, SmartCurve, The Science of Sure and associated logos are trademarks and/or registered trademarks of Hologic, Inc., and/or its subsidiaries in the United States and/or other countries. This information is intended for medical professionals in the U.S. and international markets and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative.

www.hologic.com | info@hologic.com | +32.2.711.4680

- **Curved compression surface** applies uniform compression over the entire breast.
- Shown to **improve comfort in 93% of patients** who reported moderate to severe discomfort with standard compression paddles.²
- **Processing software** takes paddle geometry into account, so there's no impact on image quality.²
- There's **no change in positioning technique or workflow** for technologists.
- Can be used with a **MammoPad® breast cushion** for additional cushioning and increased tissue capture.

Proof. Not Promises.